## Schleifscheibenwechsel

**HSS-88** (inkl. Umrüstsatz)

**& HSS-AZ/ -ZA**

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**GB** Included in delivery/ Accessories

Included in delivery	Item no.
HSS-88 Conversion kit, Dia-Quick and foil grinding disc, incl. puller and grease	168071
Spare Dia-Quick grinding disc (partially diamond coated), incl. grease	162063
Spare foil grinding disc, 1set = 4 pieces, incl. grease	162013
Spare foil grinding disc, fine, 1 set = 4 pieces, incl. grease	162023


Accessories	Item no.
Sealing grease for trimmer, 25 ml	629002
Puller for grinding disc change	164026

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## Replacement of the Dia-Quick Grinding Disc and Foil Grinding Disc

The assembly of the HSS-88 conversion kit only happens once and is only necessary when a Dia-Quick or foil grinding disc is mounted for the first time. On later occasions the used Dia-Quick and/or foil grinding disc is simply exchanged for the new one.

### WARNING:

Disconnect the power plug before starting any maintenance work.

1. Loosen the star handles to take off the front cover.
2. Loosen the central screw of the axle turning it in a left direction, whilst holding onto the grinding wheel very firmly.
3. As and if required, loosen the push rod by putting a screwdriver behind between the push rod and the grinding wheel and take it off together with the PVC wheel.
4. On the thread which is now revealed, you firmly screw on the puller which is supplied with the replacement grinding wheel and then, by screwing in the central screw, you pull the grinding wheel down off the axle.
5. Then please remove the brass bush which is still sitting in place on the stub of the shaft and exchange this for the somewhat longer brass bush from the conversion kit. For easier removal of the old brass bush it may be better to briefly push a screwdriver behind it.
6. Clean the stub of the shaft and lightly grease it with the lubricant supplied with the Dia-Quick grinding wheel / foil wheel. Now set on the support plate from the conversion kit with a light turning movement and then lay on the Dia-Quick / foil grinding wheel. Set on the push rod with PVC underlay wheel and screw really tightly.
7. If possible have the grinding wheel held in two hands and then turn the fixing screw until it is really tight.
8. Clean the contact surfaces of the casing components and seal them with the grease supplied by coating with a line of grease on the casing component and without allowing any time to elapse.
9. Put back the front plate and fasten the star handles tightly. First the bottom star handles, then the top star handles. Now spread the grease and remove any grease.

Please also use the delivered grease for greasing the windings!

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Wassermann Dental-Maschinen GmbH