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**Instructions for the assembly of the conversion kit and of the  
Dia-Quick grinding wheel / foil wheel on the  
plaster model trimmer HSS-88**

1. Loosen the turn locks to take off the front plate.
2. Loosen the central screw of the axle turning it in a left direction, whilst holding onto the grinding wheel very firmly.
3. As and if required, loosen the push rod by putting a screwdriver behind between the push rod and the grinding wheel and take it off together with the PVC wheel.
4. On the thread which is now revealed, you firmly screw on the puller which is supplied with the replacement grinding wheel and then, by screwing in the central screw, you pull the grinding wheel down off the axle.
5. Then please remove the brass bush which is still sitting in place on the stub of the shaft and exchange this for the somewhat longer brass bush from the conversion kit. For easier removal of the old brass bush it may be better to briefly push a screwdriver behind it.
6. Clean the stub of the shaft and lightly grease it with the lubricant supplied with the Dia-Quick grinding wheel / foil wheel. Now set on the support plate from the conversion kit with a light turning movement and then lay on the Dia-Quick / foil grinding wheel. Set on the push rod with PVC underlay wheel and screw really tightly.
7. If possible have the grinding wheel held in two hands and then turn the fixing screw until it is really tight.
8. Clean the contact surfaces of the casing components and seal them with the grease supplied by coating with a line of grease on the casing component and without allowing any time to elapse.
9. Put back the front plate and make sure the turn locks are really tight, removing any superfluous grease.

The assembly of the conversion kit only happens once and is only necessary when a Dia-Quick or foil grinding wheel is mounted for the first time. On later occasions the used Dia-Quick and/or foil grinding wheel is simply exchanged for the new one.

Text from following drawing:

This drawing may not be reproduced, nor be made accessible to third parties, nor may it be used by the recipient or a third party in a manner representing misuse. Infringements will lead to claims for compensation for damages and can have legal consequences.

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