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Instructions for the changing of corundum and diamond grinding wheels on
plaster model trimmer HSS-88

Pull the mains plug before changing the discs!

1. Loosen the star handles to take off the front cover.
2. Loosen the central screw on the axle counter-clockwise and hold the grinding wheel.
3. If necessary loosen the thrust piece by sticking a screwdriver between the thrust piece and the grinding wheel and take it off together with the PVC washer.
4. Screw on the puller which is supplied with the replacement grinding wheel onto the now visible thread. Then you pull off the grinding wheel down from the axle by screwing in the central screw.
5. Clean the shaft stub and brass bushing thoroughly and grease them slightly with the lubricant supplied with the replacement grinding wheel. Put on the grinding wheel by slightly turning. Put on the thrust piece and the PVC washer and screw on the central screw.
 - Check if the disc runs smooth by pushing it with the hand.
 - Advice: Use a piece of chalk and hold it steady to the disc while you twist it with the hand. Through this you could notice a possible hit which you could fix by pressing with the heel of hand.
 - Attention!: Please do not use any tool for this! In case you cannot correct the hit please contact our service team.
6. Hold the grinding wheel firm and screwing on the fixing screw tight. Check again for any hits. If you notice any hits please follow step 5.
7. Clean the contact surfaces of the housing parts and seal them with the grease supplied by applying a line of grease on one part of the housing. Do not spread it.
8. Put back the front plate and fasten the star handles tightly. Remove any, removing any spare grease.

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