SCHLEIFGERÄT KMS

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
Abb./Fig. 2

Grinding Equipment KMS

for bite orientated trimming of plaster jaw models

The production of jaw models belongs to the daily routine of a dental technician in the orthodontist field. The demand that all skull levels should be taken into consideration led automatically to the construction of an appropriate piece of equipment.

Introduction

Jaw models which are precisely prepared by the dental technician give the orthodontist the basis for the necessary therapy.

An orthodontist's diagnosis can only be given according to the anatomical facts of the upper and lower jaw as well as the intercusp and the intermaxilla relationship to one another. Also, demand for more notice to be taken of the fazilan soft parts ratio has also influenced the construction of the diagnostic models.

Because of these reasons the basic design was chosen according to „BEGG“.

Due to various angles of the upper and lower jaw molar area there are no problems when a cross over bite exists, in the area of the vestibulums.

Also, the diagnosis is easier because the dysgnathien is immediately noticeable.

The requested task, a rational and accurate method of producing plaster models is fully possible with the grinding equipment and the directional lamp which belongs to it.

All skull levels are exactly taken into consideration. Through the chosen form of the upper and lower jaw model in the vestibular area, it is possible to reach an exact assessment of the soft parts.

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
Abb./Fig. 4

contact with water are made of stainless steel. All linings and shafts are sealed with gaskets.

The working surface (3) has large slits to allow the water to flow through to avoid the models sticking during grinding.

The angle limit stop (4), moved by a thrust ball bearing, is easily pivoted away for the production of the normal plaster casts.

The limit stop for grinding the various base edges, is positioned by a precise locking device. A wedge fixes into a corresponding groove. Due to the mounted counterweight the wedge falls easily into the given angular position of the limit stop.

The installation of the plaster grinding equipment is very easy. The existing working surface is replaced by the plate of the grinding appliance (3.1).

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

Abb./Fig. 3

Abb./Fig. 3.1

Description: Grinding Equipment KMS

A simple and exact usage was the foremost thought in the construction of the equipment. It is also designed for continuous operation.

All parts that come into

Directional Lamp GTE

The directional lamp is made of stainless steel and the transformer is fully insulated (compound-filled). The halogen bulb (10 W) burns only when the plaster grinding machine is switched on. The lamp shining through an integrated optic projects a vertical ghost line onto the grinding table. Consequently, it is possible

to exactly grind the tuber level of the upper jaw by aligning the Raphe-Median-Line to the ghost line (5). In addition, the working surface in front of the trimmer grinding disc is extremely well illuminated.

The directional lamp must be installed by the manufacturer or the licensed supplier of the plaster trimmer otherwise the validity of the guarantee will be lost.

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
Abb./Fig. 7

Mould Plate for the Upper Jaw

The mould makes it easy to grind the combined tuber levels of the upper and lower jaws. The base form of the ground upper jaw fits exactly into the mould which is then set against the angle limit stop, always exactly guided towards the grinding disc.

The lower jaw laid in occlusion on the upper jaw, will thus be precisely adjusted when ground.

A continual inspection of the joint tuber levels during grinding is no longer necessary (6).

The mould can also be used for very large upper jaws when the side teeth area has to be ground (7).

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Grinding the Models

Upper Jaw

Lock the angle limit stop in the centre (90°). Lay the upper jaw model with the occlusion against the limit stop and holding it in this position (8) push it against the trimmer disc.

After grinding the base just a short while, the model pushes itself against the padded limit stop and it is no longer necessary to hold it. Remove as much as required of the base thickness.

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Mark the Raphe-Median-Line with a pencil and lay the model with the ground base onto the working surface. Move it towards the trimmer disc so that the ghost line and the Raphe-Median-Line mark correspond (10). When required dorsal limit is reached the grinding process is completed.

To grind the sides of the upper jaw turn and lock the angle limit stop in the 62° position and lay the tuber level against the limit stop (11).

Remove the plaster on both sides to the lowest point of the vestibule.

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Afterwards, symmetrically grind the upper jaw in the 25° position of the angle limit stop from the middle of the eye tooth towards the middle of the upper jaw (12).

Make the tuber level slanting by setting the angle limit stop at 90° parallel to the side teeth area (13). The upper jaw plaster model is now completed.

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
Abb./Fig. 15

Upper Jaw / Lower Jaw

Using the wax impression bring the lower jaw model in occlusion with the already ground upper jaw. Put the upper jaw model in the mould and lay it against the 90° limit stop (14).

By shifting the limit stop, grind the rear side of the lower jaw model to meet the upper jaw model. When the joint tuber level is reached the grinding noise alters considerably.

Lay the tuber level of both models, which are still in the occlusa position, on the grinding table and put the upper jaw base against the 90° limit stop. Now grind the lower jaw model down to the desired level.

The bite orientated position of the model is now reached.

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
Abb./Fig. 17

Lower Jaw

In the side position, grind the lower jaw model the same as the upper jaw. However, the angle limit stop must be in the 55° position while the base must be ground according to „BEGG“ (16).

Grind the front area using the 90° limit stop, hence parallel to the tuber level or symmetrically rounded off to the middle 3-3 (17).

Again, the last operation is to make the back edge slanting, whereby the already adjusted 90° limit stop for positioning the side teeth area is used.

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When both the upper and lower jaw models need to be set-up together e.g. for digital photography, it is recommended that the back edge angle is ground with the help of the template. Stand the template on its thin edge and position the side teeth area of the upper jaw against it (18). There is then enough space for the side teeth of the lower jaw which then lie in the template's cutout section.

- The grinding operation is finished when the bases of the upper and lower jaw models are checked.
- All angles and edges must be symmetrical to one another.
- Eventual corrections/alterations must be done by grinding.

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

Abb./Fig. 18

Features

When the above instructions are followed step-by-step, it is possible for every technician to produce exact models. After a short trial period the working time is noticeably shorter than with the usual method.

Important

The grinding machine KMS should, if possible, only be used in conjunction with a plaster trimmer that has a diamond disc. This improves the grinding process.

The technician uses less energy because the diamond disc and the plaster model trimmer do the grinding and not the technician.

Plaster Model Grinder HSS 88, HSS 88 Super, HSS 99.

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Specification and design are subject to change without notice.

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Dental-Maschinen GmbH