



MADE IN GERMANY

Wassermann is a medium-sized, owner-operated company with head office in Hamburg. We design and produce high-quality technical devices in Germany. Our main users are dental laboratories, dental clinics and various industrial areas worldwide.

Why a Wassermann product?

Wassermann products are manufactured from high-quality materials and the proportion of plastics used is intentionally kept as low as possible. Our devices are known for their suitability for daily use, functionality and longevity. As a result, the users place high levels of trust in the Wassermann brand. We feel it is our duty to support our customers on a daily basis when it comes to safety, service and productivity. Technological perfection, optimum quality and especially sustainability are the benchmarks we set. Thanks to the combination of these values, the family-owned company, since its establishment in 1927, has gained worldwide recognition as a manufacturer of high-quality products.



Further information are available on our website at www.wassermann.hamburg



/WassermannDentalmaschinen

Your Dental Dealer _____

Wapo-Ex 8 III



Boil Out Unit



Wassermann Dental-Maschinen GmbH
Rudorffweg 15-17 · 21031 Hamburg · Germany
Phone: +49 (0)40 730 926-20/ -24 · Fax: +49 (0)40 730 37 24
export@wassermann-dental.com · www.wassermann.hamburg
www.facebook.com/WassermannDentalmaschinen

Boil Out Unit **Wapo-Ex 8 III**

The **Wapo-Ex 8 III** enables the automatic boil out of up to 8 flask halves, simultaneously and fully automatically. The relevant parameters are at once visible on the touch screen, the input can be made intuitively. The user can adjust the numerous parameters and options (such as week timer) to his individual requirements by simple programming.

Splash protection
made of stainless steel

Energy saving
thermal insulation



Integrated hand spray for
manual scalding

Touch screen for
fast and easy operation

High quality construction,
most parts made of
stainless steel

Major Features

- Compact tabletop automat for scalding up to 8 flask halves
- Thermostat protected, gasketless immersion pump in bronze, maintenance-free
- Counter-sunk spray system inside of the side panels, for easy access to flask baskets
- Visual and audible water level indicator
- Easy to clean without tools
- **Included in delivery:**
Service set (4 filter mats, 1 L lime remover, 2 L wax remover)
Discharge hose, 1.5 m
Hose clip Ø 20-32 x 12.2

Accessories	Item no.
Mobile stand (W 620 x H 550 x D 430 mm)	170120
Service set (4 filter mats, 1 L lime remover, 2 L wax remover)	170021
Discharge hose per meter	320030
Wapo-Ex 8 III	Item no. 170981
Voltage	220–240 V / 50/60 Hz
Power consumption	9 A
Output	2000 W
W x H x D	610 x 540 x 420 mm
Weight	42.0 kg
Sound level	≤ 70 dB (A)
Water capacity	18.5 l