



## MADE IN GERMANY

Wassermann is a medium-sized, owner-operated company with head office in Hamburg. We design and produce high-quality technical devices in Germany. Our main users are dental laboratories, dental clinics and various industrial areas worldwide.

### Why a Wassermann product?

Wassermann products are manufactured from high-quality materials and the proportion of plastics used is intentionally kept as low as possible. Our devices are known for their suitability for daily use, functionality and longevity. As a result, the users place high levels of trust in the Wassermann brand. We feel it is our duty to support our customers on a daily basis when it comes to safety, service and productivity. Technological perfection, optimum quality and especially sustainability are the benchmarks we set. Thanks to the combination of these values, the family-owned company, since its establishment in 1927, has gained worldwide recognition as a manufacturer of high-quality products.



Further information are available on our website at [www.wassermann.hamburg](http://www.wassermann.hamburg)



/WassermannDentalmaschinen

Your Dental Dealer \_\_\_\_\_

## Wapo-Ex



### Boil Out and Polymerisation Unit



Wassermann Dental-Maschinen GmbH  
Rudorffweg 15-17 · 21031 Hamburg · Germany  
Phone: +49 (0)40 730 926-20/ -24 · Fax: +49 (0)40 730 37 24  
export@wassermann-dental.com · www.wassermann.hamburg  
www.facebook.com/WassermannDentalmaschinen

# Combination Automat Wapo-Ex

Boil Out and Polymerisation Unit

The **Wapo-Ex** is the cost-effective solution for boil out or polymerisation. The device stands out with its efficiency and reliability. Thanks to its separate compartments, a prompt switch between the two modes is possible. For quick and reproducible work processes, of course exactly to the minute and degree.

Splash protection  
made of stainless steel

Energy saving  
thermal insulation



Integrated hand spray  
for manual scalding

Informative display,  
enables a variety  
of functions  
in applications  
and in reproducible  
work sequences

High quality construction,  
most parts made  
of stainless steel

## Major Features

- Combined table automat for scalding or polymerisation
- Capacity: 4 flask halves for scalding, 6 flasks for polymerisation
- Clock timer settings accurate to the minute and to a specific degree
- Reproducible work processes via programmable scalding or polymerisation sequences
- Special function "scalding during polymerisation": Scalding is possible at the current temperature of the polymerisation process underway
- Visual and audible water level indicator
- **Included in delivery:**  
Service set (10 filter mats, 1 L lime remover, 2 L wax remover)  
Discharge hose, 1.5 m  
Hose clip Ø 20-32 x 12.2

Accessories	Item no.
Mobile stand (W 560 x H 550 x D 410 mm)	170100
Service set (10 filter mats, 1 L lime remover, 2 L wax remover)	170020
Discharge hose per meter	320030

  

Wapo-Ex	Item no. 170996
Voltage	220-240 V / 50/60 Hz
Power consumption	9 A
Output	2000 W
W x H x D	550 x 570 x 400 mm
Weight	35.5 kg
Sound level	≤ 70 dB (A)
Water capacity	approx. 20 l