

**MADE IN GERMANY**

Wassermann is a medium-sized, owner-operated company with head office in Hamburg. We design and produce high-quality technical devices in Germany. Our main users are dental laboratories, dental clinics and various industrial areas worldwide.

**Why a Wassermann product?**

Wassermann products are manufactured from high-quality materials and the proportion of plastics used is intentionally kept as low as possible. Our devices are known for their suitability for daily use, functionality and longevity. As a result, the users place high levels of trust in the Wassermann brand. We feel it is our duty to support our customers on a daily basis when it comes to safety, service and productivity. Technological perfection, optimum quality and especially sustainability are the benchmarks we set. Thanks to the combination of these values, the family-owned company, since its establishment in 1927, has gained worldwide recognition as a manufacturer of high-quality products.



Further information are available on our website at [www.wassermann.hamburg](http://www.wassermann.hamburg)



/WassermannDentalmaschinen

Your Dental Dealer \_\_\_\_\_



Wasi-Steam 2,  
basic equipment



Wasi-Clean water softener,  
available as accessory



Splash protection entirely  
made of stainless steel,  
available as accessory



Steam gun

**Wasi-Steam 2**



**Steam Cleaner**



Wassermann Dental-Maschinen GmbH  
 Rudorffweg 15-17 · 21031 Hamburg · Germany  
 Phone: +49 (0)40 730 926-20/ -24 · Fax: +49 (0)40 730 37 24  
 export@wassermann-dental.com · www.wassermann.hamburg  
 www.facebook.com/WassermannDentalmaschinen

# Steam Cleaner **Wasi-Steam 2**

The powerful **Wasi-Steam 2** is a reliable steam cleaner with a wide range of application. This steam cleaner is one of the most powerful units on the market. Thorough and efficient it is suitable for all cleaning jobs in the lab. Its high cleaning pressure of up to 6 bar also guarantees a high cleaning temperature.

Control knobs to regulate wet steam saturation and steam quantity

Clearly laid-out LED-display: display and check of the operating states

Stainless steel housing

Comfortable usage because of single-handed operation of the steam gun

Selector for wet or dry steam

**Accessory:**  
Splash protection entirely made of stainless steel

## Major Features

- User-friendly operation
- It safes gaskets because it reduces pressure to 5 bar after 10 minutes of non-use
- Permanent water supply
- Electronic delimer without chemical additives
- Service interval display
- Huge stainless steel boiler
- **Included in delivery:**  
Discharge hose  
Supply hose  
Square socket spanner

Accessories	Item no.
Wasi-Clean water softener, table top version W 190 x H 405 x D 260 mm	174995
Splash protection, stainless steel W max. 555 x H max. 250 x D max. 160 (left), 210 (right) mm	174230
Wall mount incl. fastening set	174280
Console for splash protection (for wall mounting)	174288
Special nozzle Ø 2.5 mm (Basic nozzle Ø 1.2 mm, item no. 174100)	174352
<b>Wasi-Steam 2</b>	<b>Item no. 174992</b>
Voltage	220–240 V / 50/60 Hz
Power consumption	8.3 A (Si 10 A)
Output	1880 W
W x H x D	450 x 490 x 265 mm
Weight	30 kg
Sound level	≤ 70 dB (A)
Working pressure	max. 6 bar (stand by 5 bar)
Steam temperature	max. 164°C at 6 bar
Boiler capacity	2.2 l
Safety class	IP 33
Duty	Interrupt operation 5/5 min