


● GRINDING DISCS

Well-rounded from every angle



Corundum grinding disc

Foil grinding disc

Diamond grinding disc

Dia-Quick grinding disc,
partly diamonded

Technical information

Experts know: Wassermann discs set standards. Available for every plaster hardness and offering the user top quality to meet the highest demands. From the versatile foil grinding disc to Wassermann's own diamond grinding discs – our experience is your security. Don't settle for less.

Foil grinding disc for HSS-AZ/-ZA and HSS-88

- silicon carbide coating – both sides usable
- flexible, stable design
- electronically balanced backing
- disc not glued to backing, so easily replaced
- fast and easy replacement of work discs

Corundum grinding disc for HSS-88

- approx. 10 mm thick silicon carbide coating
- self-regenerating so even, powerful cutting over a long period
- whole surface glued to backing for maximum safety
- electronically balanced to no vibration

Diamond grinding discs, see reverse side


● Model trimmers including grinding discs

HSS-99 incl. diamond grinding disc	Item no.: 164998	 HSS-99  HSS-88  HSS-AZ  HSS-ZA	Grinding disc	Item no.:
HSS-88 incl. diamond grinding disc	Item no.: 168998		Spare diamond grinding disc for HSS-99, entirely diamonded, incl. puller and grease	164063
HSS-88 incl. Dia-Quick grinding disc, partly diamonded	Item no.: 168996		Spare diamond grinding disc for HSS-88, entirely diamonded, incl. puller and grease	168063
HSS-88 incl. corundum grinding disc	Item no.: 168999		Spare Dia-Quick grinding disc for HSS-88 / -AZ / -ZA, partly diamonded, incl. grease	162063
HSS-88 incl. foil grinding disc	Item no.: 168991		Spare corundum grinding disc for HSS-88, 10 mm thick, incl. puller and grease	168024
HSS-AZ incl. Dia-Quick grinding disc, partly diamonded	Item no.: 162997		Spare foil grinding disc for HSS-88 / -AZ / -ZA, 1 set = 4 pieces, incl. grease	162013
HSS-AZ incl. foil grinding disc	Item no.: 162996		Spare foil grinding disc, fine, for HSS-88 / -AZ / -ZA, 1 set = 4 pieces, incl. grease	162023
HSS-ZA incl. Dia-Quick grinding disc, partly diamonded	Item no.: 161996	Conversion kit for HSS-88 for Dia-Quick- and foil grinding disc, incl. sealing grease and puller for grinding disc change	168071	
HSS-ZA foil grinding disc	Item no.: 161997			

Wassermann diamond-impregnated discs make easy work of even the hardest plaster – minimal pressure, maximum time saving.

Dia-Quick grinding disc for HSS-AZ/-ZA and HSS-88

- partly diamonded, island shape design
- maximum safety thanks to electronically balanced backing



Dia-Quick grinding disc, partly diamonded

Diamond grinding disc for HSS-88 and HSS-99

- full diamond impregnation for silky smooth surfaces – model surfaces are always flat despite the aggressive grinding power
- rough and fine grinding in one process
- extremely long service life – diamonds cover both sides so cannot break off
- outstanding usage option in orthodontics

Changes reserved

www.wassermann-dental.com