

● STEAM-JET WASI-STEAM CLASSIC

Minimal space requirement, maximum benefit!



Technical information

Compact format, strong performance: the Wasi-Steam classic is the perfect solution when there is little space and big demand. With its variable work pressure of up to 4 bar it offers in particular the multi-station user a gentle but quick and thorough cleaning of all small parts. The 1,200 Watt special steel heating guarantees a short heat-up phase and convincing steam temperature of up to 151°C when in use. The well designed single-hand operation of the steam gun provides maximum operating comfort.

- compact steam-jet for cleaning of all small parts
- manual filling with a volume of 3.3 litres
- optimal for use in restricted space
- comfortable application through single-handed operation of the jet gun
- variable working pressure, variable up to 4 bar
- short heat-up phase through 1,200 Watt-special steel heating
- gentle and quick to use, thorough results
- high safety through dry-run protection
- wall mounting possible

● Technical data and accessories

Wasi-Steam Classic	Item no.: 174993
Voltage*	220-240 V / 50/60 Hz
Power consumption	5.3 A (6.3 A medium-blowing-fuse)
Output	1220 W
W x H x D	235 x 340 x 340 mm
Weight	11.0 kg
Sound level	≤ 70 dB (A)
Working pressure	max. 4 bar
Steam temperature	approx. 151°C
Boiler capacity	3.3 l
Safety class	IP 22

Incl. in delivery	Item no.:
Fastening set for wall mounting	174358

Accessories	Item no.:
Jet Ø 2.5 mm	174352
Splash protection (W 235 x H 385 x D 294 mm)	174360
lime remover, 1000 ccm	655002

* Other voltages on request

Technical changes reserved



Steam gun

Splash protection

Major features at a glance

- compact steam-jet for all small parts
- easy to use, meets high safety standards
- durable, versatile and ecofriendly
- housing plastic-powder-coated

www.wassermann-dental.com