

## ● POLYMERISATION UNIT

**WAPO-MAT III**

Small size, big performance

**Technical information**

No-compromise, compact and versatile: The Wapo-Mat requires very little space, offers powerful performance and also has a practical multifunction touch-screen terminal. The water bath can be used for long-time and boiling polymerisation of all common plastics. Optimised water temperatures are guaranteed by its powerful, electronically controlled heater. Despite its compact shape, the Wapo-Mat can hold up to 9 flasks. All program steps are fully automatic, while time and temperature can be programmed to the minute/degree using the multifunctional touch display. In factory setting the main polymerization process cuts in 90 minutes automatically at 95° C 90 minutes before the program ends. The choice of 3 times and 3 temperatures is possible as an option.

- compact benchtop unit for long-time and boiling polymerisation
- fully automatic program sequence, capable to handling 9 flasks simultaneously
- multifunction touch-screen terminal with minute and degree settings, remaining time reading
- powerful, electronically controlled heater for optimal water temperature, low water level indicated visually and audibly
- energy-saving tank thermal insulation
- low maintenance, easy to clean

## ● Technical data and accessories

Wapo-Mat III	Item no.: 170980
Voltage*	220-240 V / 50/60 Hz
Power consumption	8 A
Output	1800 W
W x H x D	330 x 460 x 530 mm
Weight	18.5 kg
Sound level	≤ 70 dB (A)
Water capacity	18.5 l

Incl. in delivery	Item no.:
PVC hose 19 x 4 (3/4") – 1.5 m	320030
Hose clamp Ø 20 – 32 x 12.2	322017
Lime remover – 1000 ccm	655002

Accessories	Item no.:
PVC hose 19 x 4 (3/4") – per meter	320030
Lime remover – 1000 ccm	655002

\*Other voltages on request

Technical changes reserved



### Major features at a glance

- compact, energy-saving benchtop polymerisation unit
- fully automatic program sequence
- automatic main polymerisation
- multifunction touch-screen terminal
- simple to operate, easy to clean

[www.wassermann-dental.com](http://www.wassermann-dental.com)