

● POLISHING COMPACT UNITS WP-Ex 3000 II AND WP-Ex 2000 II

The highlight for every lab



WP-Ex 3000 II



WP-Ex 2000 II

Technical information

Powerful and durable: The WP-Ex 2000 II and WP-Ex 3000 II polishing compact units guarantee first-class polishing and extraction even in the smallest lab. Due to their compact size they can be used in the tiniest space, where they offer powerful and low-noise extraction and strong, durable and maintenance-free induction motors. The integrated daylight LEDs and the high-quality safety glass viewing protection ensure relaxed work and the right view. **Convenient:** The two soft-rubber polishing troughs could be easily removed for cleaning. The WP-Ex 3000 II has a clear multifunctional display and a multi-adjustable polishing motor. In addition the unit is equipped with a daylight LED Spot. This daylight LED spot is available as special accessory for the WP-Ex 2000 II.

- compact and powerful stainless steel polishing units for daily dental lab work
- jacket-cooled motor, dust protected
- WP-EX 2000 II: 2 settings, 1500 and 3000 min⁻¹, WP-Ex 3000 II: adjustable from 750–6000 min⁻¹
- two extra large micro fine dust filterbags, easy to replace
- impact-resistant, adjustable safety glass viewing protection
- integrated day-light LEDs
- WP-Ex 3000 II: LED spot included in delivery
WP-EX 2000 II: daylight LED spot available as special accessory
- two easily removable full rubber polishing troughs
- integrated drawer
- easy to clean

● Technical data and accessories

WP-Ex 3000 II	Item no.: 146979*
WP-Ex 2000 II	Item no.: 146978
Voltage*	220–240 V / 50/60 Hz
Power consumption	2.8 A
Output	650 W
Fan	445 m ³ /h (free-blowing)
W x H x D	680 x 486 (LED spot + 140 mm) x 527 mm* 680 x 486 x 527 mm
Weight with motor	52.0 kg
Sound level	≤ 65 dB (A)
Speed	1500 / 3000 min ⁻¹ (750–6000 min ⁻¹ adjustable in addition)*

Incl. in delivery	Item no.:
Brush spindle right	152022
Brush spindle left	152025
Supply line	592012
Micro fine dust filterbag, 1 set = 2 pieces	611005
Polishing trough, 2 pieces	2 x 146503
Drawer (W 621 x H 23 x D 219 mm)	146388
Storage bowl (W 257 x H 40 x D 198 mm)	830409
Storage bowl (W 190 x H 17 x D 150mm)	395004

Accessories	Item no.:
Brush spindle right	152022
Brush spindle left	152025
Brush spindle 125 mm right	152043
Brush spindle 125 mm left	152044
Stone holder right	152056
Stone holder left	152057
Drill chuck right 1-6 mm clamping	152049
Micro fine dust filterbag, 1 set = 2 pieces	611005
Storage bowl (W 257 x H 40 x D 198 mm)	830409
Storage bowl (W 190 x H 17 x D 150mm)	395004
Suction protective grid	146396
Daylight LED spot	146390
Fuse 6.3 A	582013

*Other voltages on request

Technical changes reserved



WP-Ex 3000 II/2000 II with a storage bowl



WP-Ex 3000 II (standard scope of delivery) equipped with a LED spot and a multifunctional display

Major features at a glance

- powerful polishing compact units
- powerful, durable and maintenance-free induction motors
- two extra large micro fine dust filterbags, easy to replace
- stainless steel housing
- low-noise
- emergency off-switch and restart protection guarantee safety at work

www.wassermann-dental.com