

● **POLISHING COMPACT UNIT****WP-Ex 10 II**

The small solution for the big demand



Technical information

Compact, powerful and really stylish: the polishing compact unit WP-Ex 10 II is the perfect solution for first class polishing and suction in the smallest space. In addition to its powerful, long-life and maintenance-free squirrel-cage motor it also convinces with a powerful and low noise suction – all included in its compact dimensions. The daylight LED illumination and the safety glass viewing protection ensure relaxed working and correct viewing angle. Very convenient for cleaning: the polishing trough made of rubber can be removed very quickly and easily. An emergency off-switch and a restart protection guarantee safety at work.

- compact and powerful polisher for daily use at laboratories even with small space
- powerful, long-life and maintenance-free squirrel-cage motor
- 2 settings, 1,500 resp. 3,000 min⁻¹
- micro fine dust filterbag, large capacity and easy to replace
- motor with jacket-cooling, dust protected
- impact-resistant, adjustable safety glass viewing protection
- integrated daylight LEDs
- daylight LED spot available as special accessory
- high quality workmanship, most parts made of stainless steel

● Technical data and accessories

WP-Ex 10 II	Item no.: 146971
Voltage*	220–240 V / 50/60 Hz
Power consumption (unit)	2.2 A
Output	500 W
Fan	445 m ³ / h (free-blowing)
W x H x D	455 x 485 x 505 mm
Weight	37 kg
Sound level	≤ 65 dB (A)
Speed	1500 / 3000 min ⁻¹

Incl. in delivery	Item no.:
Brush spindle right	152022
Polishing trough	146503
Supply line	592012
Storage bowl (W 257 x H 40 x D 198 mm)	830409
Storage bowl (W 190 x H 17 x D 150mm)	395004
Micro fine dust filterbag	---

Accessories	Item no.:
Brush spindle right	152022
Brush spindle 125 mm right	152043
Stone holder right	152056
Drill chuck right 1-6 mm clamping	152049
Storage bowl (W 257 x H 40 x D 198 mm)	830409
Storage bowl (W 190 x H 17 x D 150mm)	395004
Micro fine dust filterbag, 1 set = 2 pieces	611005
Suction protective grid, 1 pce	146396
Daylight LED spot	146390

*Other voltages on request

Technical changes reserved



WP-Ex 10 II with special accessory: daylight LED spot

Major features at a glance

- powerful polishing compact unit
- powerful, long-life and maintenance-free squirrel-cage motor
- large micro fine dust filterbag easily to replace
- low-noise extraction
- work area made of stainless steel
- easy to clean
- numerous storage options
- emergency off-switch and restart protection guarantee safety at work

www.wassermann-dental.com