

## ● SUCTION UNIT SG-1/2 D

Powerful universal suction



### Technical information

Flexible, compact and especially powerful: The SG-1/2 D is a multi-function jack-of-all-trades for broad deployment right from polishing motors up to blasting units. As is the case with most Wassermann suction units, it is designed according to the low-pressure principle and it sucks up the maximum amount of air with a very low air speed. The advantage: relaxed and virtually noiseless working conditions. A further plus point is the high suction volume whereby, even in continuous operation, the hands do not become cold. As a short-circuit rotor the ventilator motor is nearly wear-and-tear- and maintenance-free. The performance range of the SG-1/2 D is rounded off to perfection by the automatic switching-off, also of all peripheral appliances, when filter saturation point is reached, and the large volume micro fine dust filterbag.

- universal suction unit with filter saturation display
- multifunctional deployment possibilities, recommended version for the Cemat-NT blasting unit
- completely vibration-free and low-noise running
- automatic switching-off – also of hooked up appliances – at filter saturation point
- switch-over to continuous operation possible
- integrated plug for rapid power supply to connected peripheral appliances
- high-quality finish

## ● Technical data and accessories

SG-1/2 D	Item no.: 140997
Voltage*	220–240 V / 50/60 Hz
Output	105 W
Power consumption	0.5 A
Power consumption connection unit	max. 2 A
W x H x D	410 x 460 x 185 mm
Weight	12.7 kg
Sound level	61 dB(A)
Suction capacity	360 m <sup>3</sup> /h (free-blowing)

\*Other voltages on request

Incl. in delivery	Item no.:
Supply line with switchable shock proof plug	140022
Conductance sheet	140108
Micro fine dust filterbag	---
Accessories	Item no.:
Spare micro fine dust filterbag (dust class M) 1 set = 3 pieces	611015
Spare filterbag for plaster (class 1 / class 2) 1 set = 3 pieces	611017

Technical changes reserved



Fine blasting unit Cemot-NT 4  
with automatic stand  
and suction unit SG-1/2 D

### Major features at a glance

- universal suction unit according to the low pressure principle
- nearly maintenance and wear-and-tear-free
- virtually noiseless operation, no cooling down of the hands
- large volume micro fine dust filterbag
- suitable for continuous operation, high-quality finish
- comprehensive range of accessories available

[www.wassermann-dental.com](http://www.wassermann-dental.com)