

● SUCTION BOX COMPACT II

The clean solution for optimal working



Technical information

Perfect protection for first class results: The Suction Box Compact II provides a combination of maximum safety and comfort. Without the interference of dust, splinters or fumes, it allows the user to work efficiently and relaxed directly on the object. This is possible due to very large openings, a very spacious internal volume, perfect LED illumination and the special-coated protective glass screen for an optimal view. The Suction Box Compact II offers a nearly dust-free operation, and thus meets modern requirements for workplace health and safety. On request the openings for the hands can be additionally sealed with safety covers made of metal. A particular advantage are the two-sided, separate cuts for the hand-piece and the sleeves with a lot of moving area and a high standard of hygiene.

- compact and light suction box for efficient and relaxed working
- perfect protection against splinters, dust and fumes
- plenty of mobility due to large openings and spacious internal volume
- easy to open suction box with special-coated protective glass screen
- ergonomic shape for optimal physical carriage
- excellent view, first class LED illumination, high safety standard, tilt and skid resistant
- quality magnifier optionally available
- arm rests and safety covers available as accessories on request
- perfectly matched to the combination with our low-pressure suction units

● Technical data and accessories

Suction Box Compact II	Item no.: 143370	Incl. in delivery	Item no.:
Voltage*	220–240 V / 50/60 Hz	Suction Box Compact II with LED illumination and extractor nozzle in the rear	143370
Power consumption	0.1 A		
Output	25 W		
W x H x D	370 x 310 x 345 mm		
Weight	6.5 kg		
Extractor nozzle	Ø 75 mm exterior, in the rear (available on request: Ø 60 mm exterior, in the front)		
		Accessories	Item no.:
		Magnifying glass	143251
		Set of arm rests	143303
		Set of safety covers	143305
		Various adapters Ø 75 mm (inner) available on request	---

*Other voltages on request

Technical changes reserved



Suction Box Compact II with special accessories: set of arm rests



Suction Box Compact II with special accessories: set of safety covers



Suction Box Compact II

Major features at a glance

- compact and light suction box
- plenty of mobility, ideal view
- nearly maintenance-free and free of wear
- convenient accessories available
- sturdy and easy to clean, plastic powder coated

www.wassermann-dental.com