



WGS-25 Benchtop Model,
also available as Wall-mounted Model



180°-tipping mechanism



Easy and time-saving filling

MADE IN GERMANY

Wassermann is a medium-sized, owner-operated company with head office in Hamburg. We design and produce high-quality technical devices in Germany. Our main users are dental laboratories, dental clinics and various industrial areas worldwide.

Why a Wassermann product?

Wassermann products are manufactured from high-quality materials and the proportion of plastics used is intentionally kept as low as possible. Our devices are known for their suitability for daily use, functionality and longevity. As a result, the users place high levels of trust in the Wassermann brand. We feel it is our duty to support our customers on a daily basis when it comes to safety, service and productivity. Technological perfection, optimum quality and especially sustainability are the benchmarks we set. Thanks to the combination of these values, the family-owned company, since its establishment in 1927, has gained worldwide recognition as a manufacturer of high-quality products.



Further information are available on our website at
www.wassermann.hamburg

Your Dental Dealer _____

WGS-25T and WGS-25W



Plaster Storage Bins




Wassermann Dental-Maschinen GmbH
Rudorffweg 15-17 · 21031 Hamburg · Germany
Phone: +49 (0)40 730 926-0 · Fax: +49 (0)40 730 37 24
info@wassermann-dental.com · www.wassermann.hamburg

Plaster Storage Bins

WGS-25T & WGS-25W

Wall-mounted and Benchtop Modell

WGS bins keep plaster, embedding materials and pumice powder dry and protected from impurities. The electronic vibrator guarantees even distribution of all materials, without lumps forming. The established storage and dosing system is available as wall-mounted or benchtop model.



Major Features

- The established bin for all plasters, embedding materials and pumice powder
- Optimised for use with Wassermann **DSW-2** and **DSW-2/4 K** dosing scales
- **Included in delivery:**
25 kg supply container

Accessories	Item no.
25 kg supply container with lid	390030
WGS-25T, Benchtop Model WGS-25W, Wall-mounted Model	Item no. 166998 Item no. 166999
Voltage	220–240 V / 50/60 Hz
Power consumption	0.3 A
Output	70 W
W x H x D / max. D WGS-25T	445 x 1065 x 400 mm / max. 790 mm
W x H x D / max. D WGS-25W	445 x 910 x 400 mm / max. 790 mm
Weight WGS-25T	43.6 kg
Weight WGS-25W	37.8 kg
Sound level	≤ 70 dB (A)



WGS-25W, Wall-mounted Model

Other voltages on request · Technical changes reserved