



### MADE IN GERMANY

Wassermann is a medium-sized, owner-operated company with head office in Hamburg. We design and produce high-quality technical devices in Germany. Our main users are dental laboratories, dental clinics and various industrial areas worldwide.

### Why a Wassermann product?

Wassermann products are manufactured from high-quality materials and the proportion of plastics used is intentionally kept as low as possible. Our devices are known for their suitability for daily use, functionality and longevity. As a result, the users place high levels of trust in the Wassermann brand. We feel it is our duty to support our customers on a daily basis when it comes to safety, service and productivity. Technological perfection, optimum quality and especially sustainability are the benchmarks we set. Thanks to the combination of these values, the family-owned company, since its establishment in 1927, has gained worldwide recognition as a manufacturer of high-quality products.



### Wapo-Ex 8 III



### Boil Out Unit



Further information are available on our website at [www.wassermann.hamburg](http://www.wassermann.hamburg)

Your Dental Dealer \_\_\_\_\_



Wassermann Dental-Maschinen GmbH  
Rudorffweg 15-17 · 21031 Hamburg · Germany  
Phone: +49 (0)40 730 926-0 · Fax: +49 (0)40 730 37 24  
info@wassermann-dental.com · www.wassermann.hamburg

# Boil Out Unit **Wapo-Ex 8 III**

The **Wapo-Ex 8 III** enables the automatic boil out of up to 8 flask halves, simultaneously and fully automatically. The relevant parameters are at once visible on the touch screen, the input can be made intuitively. The user can adjust the numerous parameters and options (such as week timer) to his individual requirements by simple programming.

Splash protection  
made of stainless steel

Energy saving  
thermal insulation



Integrated hand spray for  
manual scalding

Touch screen for  
fast and easy operation

High quality construction,  
most parts made of  
stainless steel

## Major Features

- Compact tabletop automat for scalding up to 8 flask halves
- Thermostat protected, gasketless immersion pump in bronze, maintenance-free
- Counter-sunk spray system inside of the side panels, for easy access to flask baskets
- Visual and audible water level indicator
- Easy to clean without tools
- **Included in delivery:**  
Service set (4 filter mats, 1 L lime remover, 2 L wax remover)  
Discharge hose, 1.5 m  
Hose band clamp Ø 20-32 x 12.2

Accessories	Item no.
Service set (4 filter mats, 1 L lime remover, 2 L wax remover)	170021
Mobile stand	170120
Discharge hose per meter	320030
<b>Wapo-Ex 8 III</b>	<b>Item no. 170981</b>
Voltage	220–240 V / 50/60 Hz
Power consumption	9 A
Output	2000 W
W x H x D	610 x 540 x 420 mm
Weight	42.0 kg
Sound level	≤ 70 dB (A)
Water capacity	18.5 l