



## MADE IN GERMANY

Wassermann is a medium-sized, owner-operated company with head office in Hamburg. We design and produce high-quality technical devices in Germany. Our main users are dental laboratories, dental clinics and various industrial areas worldwide.

### Why a Wassermann product?

Wassermann products are manufactured from high-quality materials and the proportion of plastics used is intentionally kept as low as possible. Our devices are known for their suitability for daily use, functionality and longevity. As a result, the users place high levels of trust in the Wassermann brand. We feel it is our duty to support our customers on a daily basis when it comes to safety, service and productivity. Technological perfection, optimum quality and especially sustainability are the benchmarks we set. Thanks to the combination of these values, the family-owned company, since its establishment in 1927, has gained worldwide recognition as a manufacturer of high-quality products.



Further information are available on our website at [www.wassermann.hamburg](http://www.wassermann.hamburg)

Your Dental Dealer

## Wapo-Ex



### Boil Out and Polymerisation Unit



Wassermann Dental-Maschinen GmbH  
Rudorffweg 15-17 · 21031 Hamburg · Germany  
Phone: +49 (0)40 730 926-0 · Fax: +49 (0)40 730 37 24  
[info@wassermann-dental.com](mailto:info@wassermann-dental.com) · [www.wassermann.hamburg](http://www.wassermann.hamburg)

# Combination Automat Wapo-Ex-

Boil Out and Polymerisation Unit

The **Wapo-Ex** is the cost-effective solution for boil out and polymerisation. Thanks to its separate compartments, this automatic table-top vending machine allows both these processes to be executed simultaneously, and that of course exactly to the minute and degree. The device stands out with its efficiency and reliability.

Splash protection  
made of stainless steel

Energy saving  
thermal insulation

Integrated hand spray  
for manual scalding

Informative display,  
enables a variety  
of functions  
in applications  
and in reproducible  
work sequences

High quality construction,  
most parts made  
of stainless steel



## Major Features

- Combined table automat for scalding and polymerisation; simultaneous applications possible
- Capacity: 4 flask halves for scalding, 6 flasks for polymerisation
- Clock timer settings accurate to the minute and to a specific degree
- Reproducible work processes via programmable scalding and polymerisation sequences
- Visual and audible water level indicator
- **Included in delivery:**  
Service set (10 filter mats, 1 L lime remover, 2 L wax remover)  
Discharge hose, 1.5 m  
Hose band clamp Ø 20-32 x 12.2

Accessories	Item no.
Service set (10 filter mats, 1 L lime remover, 2 L wax remover)	170020
Mobile stand	170100
Discharge hose per meter	320030

Wapo-Ex	Item no. 170996
Voltage	220-240 V / 50/60 Hz
Power consumption	9 A
Output	2000 W
W x H x D	550 x 570 x 400 mm
Weight	35.5 kg
Sound level	≤ 70 dB (A)
Water capacity	approx. 20 l