



MADE IN GERMANY

Wassermann is a medium-sized, owner-operated company with head office in Hamburg. We design and produce high-quality technical devices in Germany. Our main users are dental laboratories, dental clinics and various industrial areas worldwide.

Why a Wassermann product?

Wassermann products are manufactured from high-quality materials and the proportion of plastics used is intentionally kept as low as possible. Our devices are known for their suitability for daily use, functionality and longevity. As a result, the users place high levels of trust in the Wassermann brand. We feel it is our duty to support our customers on a daily basis when it comes to safety, service and productivity. Technological perfection, optimum quality and especially sustainability are the benchmarks we set. Thanks to the combination of these values, the family-owned company, since its establishment in 1927, has gained worldwide recognition as a manufacturer of high-quality products.



Polyquick



Polymerisation Unit



Further information are available on our website at www.wassermann.hamburg

Your Dental Dealer _____



Wassermann Dental-Maschinen GmbH
Rudorffweg 15-17 · 21031 Hamburg · Germany
Phone: +49 (0)40 730 926-0 · Fax: +49 (0)40 730 37 24
info@wassermann-dental.com · www.wassermann.hamburg

Polymerisation Unit

Polyquick

The versatile **Polyquick** polymerisation unit allows manual and timer-controlled use. This automat is designed for polymerisation of up to 3 flasks in flask holder or 2 cast resin flasks or rebasing frames or 6 orthodontic regulatings at one time. Polymerisation times can be set up to 59 hours and 40 minutes. The desired polymerisation temperature can be set up to 95 °C in steps of 5. The heating is cut off automatically when the temperature is reached.



Energy saving thermal insulation

Durable and robust metal housing, plastic powder coated

Informative display with multiple functions and application options

Dry-run protection

Needle valve

Major Features

- Pressure polymerisation unit for all applications
- An interlock switch ensures that compressed air is only released when the door is completely closed
- Pressure relief valve for protection in case of malfunction (7 bar)
- Fast emptying of hot water using compressed air and needle valve
- Easy operation, low maintenance
- **Incl. in delivery:**
 - Flask bracket support
 - Flask basket
 - Compressed air hose, 2 m
 - Discharge hose, 2 m

Accessories	Item no.
Compressed air hose per meter	320026
Discharge hose per meter	320038

Polyquick	Item no. 170999
Voltage	220-240 V / 50/60 Hz
Power consumption	4.2 A
Output	1000 W
W x H x D	340 x 435 x 355 mm
Weight	21.1 kg
Sound level	≤ 70 dB (A)
Compressed air connection	min. 4 bar
Operating pressure	max. 6 bar
Pressure chamber	Ø 150 mm, depth 285 mm