



Additional tank



Gobi-3



Gobi-4

MADE IN GERMANY

Wassermann is a medium-sized, owner-operated company with head office in Hamburg. We design and produce high-quality technical devices in Germany. Our main users are dental laboratories, dental clinics and various industrial areas worldwide.

Why a Wassermann product?

Wassermann products are manufactured from high-quality materials and the proportion of plastics used is intentionally kept as low as possible. Our devices are known for their suitability for daily use, functionality and longevity. As a result, the users place high levels of trust in the Wassermann brand. We feel it is our duty to support our customers on a daily basis when it comes to safety, service and productivity. Technological perfection, optimum quality and especially sustainability are the benchmarks we set. Thanks to the combination of these values, the family-owned company, since its establishment in 1927, has gained worldwide recognition as a manufacturer of high-quality products.



Further information are available on our website at www.wassermann.hamburg

Your Dental Dealer _____

Gobi-2, Gobi-3, Gobi-4



Fine Blasting Units



Wassermann Dental-Maschinen GmbH
Rudorffweg 15-17 · 21031 Hamburg · Germany
Phone: +49 (0)40 730 926-0 · Fax: +49 (0)40 730 37 24
info@wassermann-dental.com · www.wassermann.hamburg

Fine Blasting Units

Gobi-2, Gobi-3 & Gobi-4

The **Gobi-2** enables you to work not only very efficiently and quick but is also extremely flexible. Due to its modular structure the basic unit can be extended by 1 or 2 blasting material tanks, up to 4 chambers. A unit that satisfied of long service life and highest flexibility: Sand-blasting, roughening or polishing using all common blasting materials and grain sizes.

We recommend the use of Wassermann blasting materials and our low-pressure suction units **SG-10**, **SG-1/1** and **SG-1/2 D**.

Precise working pressure of 0.8–6 bar with pressure gauge control

Modular structure, extendable to up to 4 chambers

Easy choice of the desired blasting material by using the selection switch for colour circuit settings

High quality construction, made of metal, designed for long durability

Rapid ventilation of the blasting material containers for rapid re-filling and changing of blasting material

Optimum workspace illumination, protected fluorescent lamp

Comfortably to handle with ergonomically shaped handpieces



Major Features

- Genuine quick Wassermann start/stop function, foot-controlled
- Cost savings thanks to noticeable low consumption of blasting material
- High hygiene standard because of changeable sleeves
- Special particularly durable hoses
- Convenient colour coding system
- Specially coated safety screen for ideal visibility
- With plug socket, e.g. for a suction unit
- **Included in delivery:**
Compressed air hose, 2 m
Supply line
Gobi-2 incl. 2 fine blasting nozzles: 0.8 (blue), 1.2 (red)
Gobi-3 incl. 3 fine blasting nozzles: 0.8 (blue), 1.2 (red) and optionally additional tank with fine blasting nozzle left 0.6 (yellow) or right 0.8 (green)
Gobi-4 incl. 4 fine blasting nozzles: 0.8 (blue), 1.2 (red) and additional tank, left with fine blasting nozzle 0.6 (yellow) and additional tank, right with fine blasting nozzle 0.8 (green)

Accessories

Accessories	Item no.
Additional tank, left with fine blasting nozzle 0.6 mm (yellow)	190800
Additional tank, right with fine blasting nozzle 0.8 mm (green)	190802
Compressed air hose per meter	320026
Spare fine blasting nozzle 0.6 mm	190550
Spare fine blasting nozzle 0.8 mm	190551
Spare fine blasting nozzle 1.2 mm	190552

Gobi-2

Gobi-2	Item no. 190896
Voltage	220–240 V / 50/60 Hz
Power consumption	0.2 A
Plug socket power consumption	max. 8 A
Output	50 W
W x H x D	430 x 350 x 530 mm
Weight	17.5 kg
Sound level	≤ 70 dB (A)
Tank capacity	800 cm ³ = 1.2 kg
Compressed air connection	3–6 bar
Working pressure	0.8–6 bar
Extractor nozzle	Ø 75 mm exterior