



MADE IN GERMANY

Wassermann is a medium-sized, owner-operated company with head office in Hamburg. We design and produce high-quality technical devices in Germany. Our main users are dental laboratories, dental clinics and various industrial areas worldwide.

Why a Wassermann product?

Wassermann products are manufactured from high-quality materials and the proportion of plastics used is intentionally kept as low as possible. Our devices are known for their suitability for daily use, functionality and longevity. As a result, the users place high levels of trust in the Wassermann brand. We feel it is our duty to support our customers on a daily basis when it comes to safety, service and productivity. Technological perfection, optimum quality and especially sustainability are the benchmarks we set. Thanks to the combination of these values, the family-owned company, since its establishment in 1927, has gained worldwide recognition as a manufacturer of high-quality products.



Further information are available on our website at www.wassermann.hamburg

Your Dental Dealer

Puna



Fine Blasting Units



Wassermann Dental-Maschinen GmbH
Rudorffweg 15-17 · 21031 Hamburg · Germany
Phone: +49 (0)40 730 926-0 · Fax: +49 (0)40 730 37 24
info@wassermann-dental.com · www.wassermann.hamburg

Fine Blasting Units **Puna**

Sand-blasting, roughening or polishing in the premium compact class! The Puna fine blasting unit with two tanks enables you to work extremely efficiently and quickly, reliable for many years. It can be used with all common blasting materials and grain sizes for surface finishing. The specialist also impresses thanks to its low consumption of blasting material - and it is not a matter of course in its class.

We recommend the use of Wassermann blasting materials and our low-pressure suction units **SG-10**, **SG-1/1** and **SG-1/2 D**.

Precise working pressure of 0.8–6 bar with pressure gauge control

Easy choice of the desired blasting material by using the selection switch for colour circuit settings

With plug socket, e.g. for a suction unit



Blasting material tanks made of metal, high quality construction

Optimum workspace illumination, protected fluorescent lamp

Comfortably to handle with ergonomically shaped handpieces

Major Features

- Cost-effective fine blasting unit with two tanks
- Easy operation by foot pedal
- Noticeable low consumption of blasting material due to the Wassermann air mix system
- High hygiene standard because of changeable sleeves
- Special particularly durable hoses
- Convenient colour coding system
- Specially coated safety screen for ideal visibility
- Designed for long durability
- **Included in delivery:**
Puna incl. 2 fine blasting nozzles: 0.8 (blue), 1.2 (red)
Compressed air hose, 2 m
Supply line

Accessories	Item no.
Compressed air hose per meter	320026
Spare fine blasting nozzle 0.6 mm	190550
Spare fine blasting nozzle 0.8 mm	190551
Spare fine blasting nozzle 1.2 mm	190552

Puna	Item no. 190897
Voltage	220–240 V / 50/60 Hz
Power consumption	0.2 A
Plug socket power consumption	max. 8 A
Output	50 W
W x H x D	400 x 335 x 455 mm
Weight	10.5 kg
Sound level	≤ 70 dB (A)
Tank capacity	420 cm ³
Compressed air connection	3–6 bar
Working pressure	0.8–6 bar
Extractor nozzle	Ø 75 mm exterior