



Automatic starting device AS-100, available as accessory

MADE IN GERMANY

Wassermann is a medium-sized, owner-operated company with head office in Hamburg. We design and produce high-quality technical devices in Germany. Our main users are dental laboratories, dental clinics and various industrial areas worldwide.

Why a Wassermann product?

Wassermann products are manufactured from high-quality materials and the proportion of plastics used is intentionally kept as low as possible. Our devices are known for their suitability for daily use, functionality and longevity. As a result, the users place high levels of trust in the Wassermann brand. We feel it is our duty to support our customers on a daily basis when it comes to safety, service and productivity. Technological perfection, optimum quality and especially sustainability are the benchmarks we set. Thanks to the combination of these values, the family-owned company, since its establishment in 1927, has gained worldwide recognition as a manufacturer of high-quality products.



Further information are available on our website at www.wassermann.hamburg

Your Dental Dealer

SG-1/1



Suction Unit



Wassermann Dental-Maschinen GmbH
Rudorffweg 15-17 · 21031 Hamburg · Germany
Phone: +49 (0)40 730 926-0 · Fax: +49 (0)40 730 37 24
info@wassermann-dental.com · www.wassermann.hamburg

Suction Unit SG-1/1

The 1-place-suction **SG-1/1** works with Wassermann's reliable low pressure principle; i.e. with suction at very low air speed. The advantage: Relaxed to work with because of its low noise. A further advantage is the high suction volume, where the hands will not become cold even in continuous use. The fan motor is a maintenance-free squirrel-cage motor, suitable for continuous operation and almost free from wear.

Major Features

- 1-place-suction that works with the low pressure principle, multifunctional use – from polishing trough to fine blasting units
- Low noise operation, suitable for continuous operation
- Wide range of accessories available, recommended special accessory SG-1/1: Automatic starting device **AS-100** for the automatic start of suction units
- **Included in delivery:**
Supply line with switchable shock proof plug
Micro fine dust filterbag (yellow), dust class M

Accessories	Item no.
Automatic starting device AS-100	143100
Spare micro fine dust filterbag (yellow), dust class M, 1 set = 4 pieces	611010
Spare filterbag for plaster (white), class 1 / class 2 1 set = 4 pieces	611018
SG-1/1	Item no. 140995
Voltage	220–240 V / 50/60 Hz
Power consumption	0.75 A
Plug socket power consumption	max. 8 A
Output	170 W
W x H x D	260 x 400 x 330 mm
Weight	9.9 kg
Sound level	59 dB(A)
Suction capacity	500 m ³ /h (free-blowing)
Extractor nozzle	Ø 75 mm exterior

Safety and health in the workplace:
The SG-1/1 is equipped with a high-quality, easily to replace filterbag which is optionally available for micro fine dust or for plaster

Maximum flexibility because of compact construction design

Long-life brushless motor

With plug socket, e.g. for laboratory technology machines or polishing motors

