

● **WAX SCALDING****AND POLYMERISATION UNIT WAPO-EX**

The perfect everyday combination

**Technical information**

A strong duo: The Wapo-Ex from Wassermann is an all-in-one answer to your wax-scalding and polymerisation units. What is thereby so attractive about the Wapo-Ex is in particular its efficiency. Thanks to separate chambers it is possible to use this table unit for parallel scalding and hardening off – timed down to the exact minute and/or degree. The sophisticated multifunction display enables a variety of functions to be carried out both in applications and in reproducible work sequences by means of its simple-to-use programming. In addition the Wapo-Ex has optimal thermo-insulation of the containers and is consequently especially energy-saving. And to complete things, the stainless steel finish guarantees durability and low levels of wear-and-tear.

- combined table unit for scalding and polymerisation; simultaneous applications possible
- flexible and informative multi-function display with clock timer settings accurate to the minute and to a specific degree
- capacity: 4 flask halves for scalding, 6 flasks for polymerisation
- in factory setting the main polymerization process cuts in 90 minutes automatically at 95° C 90 minutes before the program ends - the choice of 3 times and 3 temperatures is possible as an option.
- reproducible work sequences via programmable scalding sequences and a standard programme
- integrated adjustable showerhead for manual scalding
- energy-saving thermo-insulation
- optical and acoustic low-water display

● Technical data and accessories

Wapo-Ex	Item no.: 170996
Voltage*	220-240 V / 50/60 Hz
Power consumption	9 A
Output	2000 W
W x H x D	550 x 570 x 400 mm
Weight	35.5 kg
Sound level	≤ 70 dB (A)
Water capacity	approx. 20 l

Incl. in delivery	Item no.:
Service set (filters, lime remover, wax remover)	170020
PVC hose 19 x 4 (3/4") – 1.5 m	320030
Hose band clamp Ø 20 - 32 x 12.2	322017

Accessories	Item no.:
Service set (filters, lime remover, wax remover)	170020
Mobile base for Wapo-Ex	170100
PVC hose 19 x 4 (3/4") – per meter run	320030

*Other voltages on request

Technical changes reserved



Major features at a glance

- combined unit for scalding and polymerisation
- high efficiency via the possibility of parallel application
- multi-purpose function and application possibilities
- simple to handle
- remarkably durable
- casing made of stainless steel, easy maintenance

www.wassermann-dental.com