

● Power vibrators KV-16 / KV-26 KV-26 plus / KV-36 and KV-56

Well vibrated means well compacted



KV-16, KV-26 plus, KV-36
stainless steel version plastic-powder-coated



KV-26
stainless steel version plastic-powder-coated

Technical information

These are real winners – solid, strong and sturdy. Powerful vibrators from Wassermann can meet any challenge. Their centrally located oscillating magnets and high manufacturing accuracy guarantee even vibration over the entire table as well as powerful vibration intensity. All machines are designed for continuous operation and can be switched from 3000 to 6000 oscillations – some with infinitely adjustable vibration intensity. All units are equipped with vibration damping rubber feet.

- extremely powerful vibrator for even and void-free shaking
- powerful vibration intensity thanks to the highest manufacturing accuracy
- switchable from 3000 to 6000 oscillations, and infinitely adjustable vibration intensity, depending on the model
- even vibration intensity across the entire table, thanks to centrally mounted magnets suitable for continuous loads
- removable and washable rubber supports
- all machines have a stainless steel yoke on the left-hand side which can be used when removing casts

● Technical data and accessories

KV	stainless steel white	stainless steel matt
KV-16	Item no.: 169998	169986
KV-26	Item no.: 169996*	169987*
KV-26 plus	Item no.: 169992**	169988**
KV-36	Item no.: 169999***	169989***
Voltage*	220-240 V / 50/60 Hz	
Power consumption	0.3 A, 0.6 A***	
Output	70 W, 130 W***	
W x H x D	145 x 150 x 143 mm 228 x 150 x 143 mm*/** 308 x 173 x 208 mm***	
Weight	3.9 kg, 5.2 kg*/**, 10.9 kg***	
Sound level	≤ 59 dB (A)	
Vibrating intensity	3000 / 6000 vibr./min (infinitely variable **/***)	
Table dimensions	130 x 130 mm, 210 x 130 mm*/** 293 x 200 mm***	
Safety class	IP 54	

*Other voltages on request

KV-56	Item no.: 169997
Voltage*	220-240 V / 50/60 Hz
Power consumption	1.3 A
Output	282 W
W x H x D	445 x 187 x 286 mm
Weight	18.6 kg
Sound level	≤ 70 dB (A)
Vibrating intensity	3000 / 6000 vibr./min. (infinitely variable)
Table dimensions	430 x 286 mm
Safety class	IP 54

Accessories	Item no.:
Spare rubber plate KV-16	169201
Spare rubber plate KV-26 / KV-26 plus	169101
Spare rubber plate KV-36	169010
Spare rubber plate KV-56	169510

Technical changes reserved



KV-16, KV-26 plus, KV-36 stainless steel version matt



KV-26 stainless steel version matt



KV-56 plastic-powder-coated cast housing

Major features at a glance

- powerful vibrator in different versions
- KV-16, KV-26, KV-26 plus and KV-36: extremely robust housing in stainless steel plastic-powder-coated or matt
- KV-56: plastic-powder-coated cast housing
- heavy weight for maximum stability
- for even and void-free shaking
- some with infinitely adjustable vibration intensity
- designed for continuous loads
- all units with vibration damping rubber feet

www.wassermann-dental.com