

“developed
for more
flexibility”

- **efficient & powerful**

Our Proven Bestseller:

Wassermann's best-selling trimmer
- now even stronger, even more
efficient, even more flexible.

OUR PROVEN RECIPE FOR SUCCESS:

We combine technology, experience,
competence and team spirit.



Always available

Would you like to request a quote or
place an order? Or do you have
questions about your order, delivery
or warranty?

We look forward to your call and
will be happy to help you.



www.wassermann.hamburg



fb.com/WassermannDentalmaschinen

**Dry Trimmer
HSS-TR**



Service Hotline

+49 (0)40 730 926-0



SCAN ME



Dry Trimmer HSS-TR

NEW Powerful
in a new look


QUALITY MADE IN GERMANY SINCE 1927

Its unique, central connection for the suction system guarantees the best suction performance. Due to its powerful, high-torque but also energy-efficient motor highest grinding speeds are achieved.



The most functions and features at a glance:

- Versatile: Various Wassermanns grinding wheels can be used
- Very high-torque and maintenance-free motor (at max. 3000 rpm)
- Stable and solid stainless steel and noble aluminum
- Central connection for suction system
- Powerful disc suction and internal suction
- Safe and easy to open, no unnecessary seals
- Large work area of 10,5 x 13,5 cm
- Grinding jigs KMS are also available in Swiss and EOS type.

Notes:

Technical Data:

Output: 1.000 W
Power consumption: 4,6 A

Weight without grinding disc: 19,8 kg
W x H x D 330 x 335 x 410 mm

Speed: 3000 U/min
Voltage: 220–240 V

Support Service:

Wassermann Dental-Maschinen GmbH
Rudorffweg 15 - 17
21031 Hamburg

Tel: +49 (0)40 730 926-0
Fax: +49 (0)40 730 37 24
www.wassermann-dental.com
sales@wassermann-dental.com