



Stainless steel clamp clips



Universal flask clamps



WW-33,  
version for 1-3 flasks

**MADE IN GERMANY**

Wassermann is a medium-sized, owner-operated company with head office in Hamburg. We design and produce high-quality technical devices in Germany. Our main users are dental laboratories, dental clinics and various industrial areas worldwide.

**Why a Wassermann product?**

Wassermann products are manufactured from high-quality materials and the proportion of plastics used is intentionally kept as low as possible. Our devices are known for their suitability for daily use, functionality and longevity. As a result, the users place high levels of trust in the Wassermann brand. We feel it is our duty to support our customers on a daily basis when it comes to safety, service and productivity. Technological perfection, optimum quality and especially sustainability are the benchmarks we set. Thanks to the combination of these values, the family-owned company, since its establishment in 1927, has gained worldwide recognition as a manufacturer of high-quality products.



Further information are available on our website at [www.wassermann.hamburg](http://www.wassermann.hamburg)



/WassermannDentalmaschinen

Your Dental Dealer \_\_\_\_\_

**WW-33**



**Hydraulic Press**



Wassermann Dental-Maschinen GmbH  
Rudorffweg 15-17 · 21031 Hamburg · Germany  
Phone: +49 (0)40 730 926-20/ -24 · Fax: +49 (0)40 730 37 24  
export@wassermann-dental.com · www.wassermann.hamburg  
www.facebook.com/WassermannDentalmaschinen

# Hydraulic Press WW-33

The WW-33 is designed for daily use. It has a working pressure that can be adjusted up to 6 tonnes – for the appropriate pressing of 1–3 or 2–4 flasks. It works efficiently thanks to the valve which can be opened to allow the piston to return to its starting position automatically – hence the press is ready again for operation without losing time.



Rust-proof piston  
made of super-strong material

Pressure infinitely adjustable  
up to 6 tonnes

Quality pressure gauge  
with precision display  
recessed in housing

Automatic piston return  
after opening of valve

Integrated  
oil level display

## Major Features

- Break-proof hydraulic press in drawn steel for pressing 1–3 or 2–4 flasks
- Distance of columns of 152 mm to allow use of larger flasks and flask clamps
- Plastic powder coated housing
- High stability
- **Included in delivery:**  
WW-33 for 2-4 or 1-3 flasks  
Pump lever

Accessories	Item no.
Universal flask	129990
Stainless steel clamp clip for 3 flasks	129992
Stainless steel clamp clip for 2 flasks	129993
Stainless steel clamp clip for 1 flask	129994
Universal flask clamp for 3 flasks	129995
Universal flask clamp for 2 flasks	129996
Universal flask clamp for 1 flask	129997
Hydraulic oil 500 ccm	629005

WW-33, 2–4 flasks WW-33, 1–3 flasks	Item no. 126998 Item no. 126999
W x D / D with pump lever	200 x 158 / 295 mm
Height WW-33, 2–4 flasks	max. 710 / min. 540 mm
Height WW-33, 1–3 flasks	max. 650 / min. 480 mm
Distance between columns	152 mm
Distance plunger plate WW-33, 2–4 flasks	max. 272 / min. 102 mm
Distance plunger plate WW-33, 1–3 flasks	max. 212 / min. 42 mm
Weight WW-33, 2–4 flasks	16.0 kg
Weight WW-33, 1–3 flasks	14.7 kg
Working pressure	up to 6 tons
Working table	Ø 95 mm