

“Engineered for Maximum Flexibility

Powerful • Flexible • Durable

Our new **Wamix Economy** is a user-friendly vacuum mixing device suitable for all types of plaster, investment materials, and silicones.

OUR PROVEN FORMULA FOR SUCCESS:

We combine technology, experience, expertise, and team spirit.



Always here for

Do you want to request a quote, order products, or have questions about your order, delivery, or warranty?

We look forward to your call and are happy to assist you.

Wamix Economy Vacuum Mixer



Info-Hotline

+49 (0)40 730 926-0



www.wassermann.hamburg



fb.com/WassermannDentalmaschinen

SCAN ME



Wamix Economy Vacuum Mixer

NEW Now easier and
faster to use

QUALITY MADE IN GERMANY SINCE 1927

The mixing device delivers excellent, bubble-free mixing results. Its simple operation and extreme durability provide reliable ease of use. A newly developed control system enables fast and straightforward operation for the user.



Key Features at a Glance:

- Designed for long-lasting durability
- Versatile as both a wall-mounted and tabletop device (table stand not included)
- Low maintenance and easy to clean
- Suitable for all dental plaster types, investment materials, and silicones
- Digital timer
- Includes mixing cup and agitator
- Angled motor for easier positioning of the mixing cup
- Convenient one-handed operation for optimal mixing results
- High-quality construction, made of metal
- Easy to assemble with just a few simple steps

Wamix Economy
Item No. 171972

Technical Specifications

Motor Power: 94 W
Current: 0.44 A

Dimensions (W x H x D): 195 x 310 x 180 mm
Vacuum: Approx. 120 mbar (absolute)

Weight: 9.3 kg
Protection Class: IP 22

Info & Consultation:

Wassermann Dental-Maschinen GmbH
Rudorffweg 15 - 17
21031 Hamburg, Germany

Tel: +49 (0)40 730 926-0
www.wassermann-dental.com
sales@wassermann-dental.com