



Recommended Retail Prices 2026



PRICE LIST 2026

This price list is effective from 01.01.2026, ex works Hamburg and is subject to change without prior notice and cannot be considered binding to the seller for any specific period.

Prices are quoted in EUR (EURO), without V.A.T. and subject to our normal terms and conditions of sale regarding delivery and payment. You could request these terms by phone or download them from our homepage.

Specifications are subject to alteration without prior notice.

Orders placed below EUR 50.-- are subject to a handling fee of EUR 10.--.

The plastic coating usually is:
natural white close to RAL 9016 and black close to RAL 9005.

Colours are also subject to change without notice. Alternative colours are available at additional surcharge.

All units are available in 220-240 V / 50/60 Hz.
Other voltages upon enquiry. Please let us know.
Units equipped with heating are only available in 220-240 V / 50/60 Hz.

Wassermann Dental-Maschinen GmbH
Rudorffweg 15 – 17, 21031 Hamburg, Germany

Export Department:
Phone: +49 40 730 926 22 or +49 40 730 926 24

export@wassermann-dental.com

www.wassermann.hamburg

www.facebook.com/WassermannDentalmaschinen

Contents

Modell production.....	3
Plaster Dispenser and Accessories.....	3
Dosing scales.....	4
Accessories for DSW-2/ -2/4K.....	4
Vacuum mixers.....	4
Accessories for vacuum mixers.....	5
Vibrators.....	5
Accessories for vibrators.....	6
Trimmers.....	6
Accessories for trimmers.....	7
Spare grinding discs.....	8
Modelling.....	9
Laboratory technology machines.....	9
Accessories for laboratory technology machines.....	10
Steam cleaners.....	10
Accessories Wasi-Steam Classic/ Classic II.....	12
Splash protectors.....	12
Boil-Out and polymerisation units.....	12
Accessories for boil out and polymerisation units.....	13
Compactors.....	14
Mixing and dosing unit.....	14
Accessories for mixing and dosing unit.....	14
Ejection press.....	15
Flask presses.....	15
Accessories for flask presses.....	16
Surface finish.....	16
Fine blasting units.....	16
Accessories for Gobi-2 (Item no. 190896).....	19
Accessories for Puna (Item no. 190897).....	19
Polishing compact units.....	19
Polishing motor.....	20
Polishing troughs.....	20
Accessories for polishing lathes.....	20
Filter for polishing compact units.....	21
Automatic stand.....	21
Suction technology.....	21
Suction Box.....	23
Accessories for suction units.....	23
Filters for suction units.....	24
Consumable materials.....	25
Duplicating silicone.....	25
Blasting materials.....	26

Price in €

Modell production**Plaster Dispenser and Accessories**

166997 D	Plaster Dispenser WGS-25W II	4.660,00
	Wall Model, rotatable (180 °), for 25 kg plaster container	
166115 D	Table-Top Support for WGS-25 W II	670,00

Dosing scales

172988 A	DSW-2 (with 1 pump, not cooled)	2.505,00
	For exact dosage of plaster and investment material	
	You can order up to 5 additional pumps	

Accessories for DSW-2/ -2/4K

172299 C	Additional pump equipment (not cooled)	229,00
-----------------	---	--------

Price in €**Vacuum mixers**

171970 A	Wamix-Touch with vacuum pump, incl. wall holder and mixing cup 350 ml with touch screen for fast and easy operation	2.245,00
171971 A	Wamix-Touch Injector with injector technique, incl. wall holder and mixing cup 350 ml with touch screen for fast and easy operation	2.030,00
171972 B	Wamix-Economy with vacuum pump, incl. wall holder and mixing cup 350 ml	1.465,00

Accessories for vacuum mixers**for Wamix-Touch/ Wamix-Touch Injector and
Wamix-2/ -2M/ -2MS/ -3M/ -2000/ CCI¹**

CCI¹ = Wamix-Comfort Plus, -Comfort, -Classic, -Injector

171295 C	Mixing cup Wamix mini, 75 ml (75 ml mixing cup not suitable for beaker mount Wamix-Comfort Plus)	219,00
171230 C	Mixing cup Wamix small, 200 ml	229,00
171283 C	Mixing cup Wamix medium, 350 ml	236,00
171237 C	Mixing cup Wamix large, 600 ml	244,00
171269 C	Mixing cup Wamix extra large, 1000 ml (1000 ml mixing cup not suitable for beaker mount Wamix-Comfort Plus)	311,00
171744 C	Table stand Wamix-Touch/ Wamix-Touch Injector	357,00
171122 C	Set of sealings (set = 5 pcs.)	9,00
610025 C	Sinterfilter	24,50

Price in €**Vibrators**

169997 A	KV-56 , <u>plastic-powder-coated</u> cast housing White Table dimensions 430 x 286 mm, vibration intensity infinitely variable	1.680,00
169999 A	KV-36 , <u>plastic-powder-coated</u> stainless steel White Table dimensions 293 x 200 mm, vibration intensity infinitely variable	810,00
169992 A	KV-26 plus , <u>plastic-powder-coated</u> stainless steel White Table dimensions 210 x 130 mm, vibration intensity infinitely variable	702,00
169996 A	KV-26 , <u>plastic-powder-coated</u> stainless steel White Table dimensions 210 x 130 mm	576,00
169998 A	KV-16 , <u>plastic-powder-coated</u> stainless steel White Table dimensions 130 x 130 mm	499,00
169989 A	KV-36 , <u>matt</u> stainless steel Table dimensions 293 x 200 mm, vibration intensity infinitely variable	810,00
169988 A	KV-26 plus , <u>matt</u> stainless steel Table dimensions 210 x 130 mm, vibration intensity infinitely variable	702,00
169987 A	KV-26 , <u>matt</u> stainless steel Table dimensions 210 x 130 mm	576,00
169986 A	KV-16 , <u>matt</u> stainless steel Table dimensions 130 x 130 mm	499,00

Accessories for vibrators

169510 C	Spare rubber plate KV-56	139,00
169010 C	Spare rubber plate KV-36	65,50
169101 C	Spare rubber plate KV-26/ KV-26 plus	56,50
169201 C	Spare rubber plate KV-16	32,50

Trimmers

164998 A	HSS-99 - 1200 Watt – incl. diamond grinding disc incl. AquaStop, large grinding area	4.975,00
168988 A	HSS-88 Super 1.500 Watt incl. Diamond grinding disc incl. AquaStop, grinding disc entirely diamonded	4.165,00
168998 A	HSS-88 – 1.000 Watt - incl. diamond grinding disc incl. AquaStop, grinding disc entirely diamonded	3.915,00
168996 A	HSS-88 – 1.000 Watt - incl. Dia-Quick grinding disc incl. AquaStop, with partially diamond disc	3.450,00
163990 B	HSS-TR Dry Trimmer Dia-Quick with partially diamond disc. Suction units on Page 22!	3.295,00
168989 A	HSS-88 Super 1.500 Watt incl. corundum-grinding disc incl. AquaStop, corundum-grinding disc	2.650,00
168999 A	HSS-88 - 1000 Watt – incl. corundum-grinding disc incl. AquaStop, corundum-grinding disc	2.390,00
168991 A	HSS-88 - 1000 Watt - incl. foil grinding disc incl. AquaStop, foil disc double-sided-use	2.390,00
162997 A	HSS-AZ - 750 Watt – incl. Dia-Quick grinding disc incl. AquaStop, with partially diamond disc	2.665,00
162996 A	HSS-AZ - 750 Watt - incl. foil grinding disc incl. AquaStop, foil disc double-sided-use	1.670,00
161996 A	HSS-ZA incl. Dia-Quick grinding disc and brush spindle on left side combination trimmer/polishing motor -500 Watt-	3.735,00
161997 A	HSS-ZA incl. foil grinding disc, double-sided-use and brush spindle on left side combination trimmer/polishing motor -500 Watt-	2.605,00
143123 C	Polishing trough without extractor nozzle incl. adjustable protection sheet and removable rubber trough	271,00

Accessories for trimmers
**Grinding attachment for HSS-99
 for producing 3D orthodontic models:**

164130 D	Model table KMS according to Begg	2.275,00
168088 D	Laser light guide for model table	928,00
164034 D	Orthodontic template, lower jaw/ upper jaw/ tuber plane	55,00

! **The Grinding attachment KMS can be delivered also as Swiss and EOS-versions. Price on demand.** **!**

164161 C	Tray for trimmer HSS-99 fitting in a sink (We always need dimensions of the sink and design 1:1)	335,00
168160 C	Tray for trimmer HSS-88/ AZ fitting in a sink (We always need dimension of the sink and design 1:1)	290,00
164037 C	Inclined bearing table for HSS-99	107,00
168037 C	Inclined bearing table for HSS-88/HSS-AZ/HSS-ZA	105,00
168168 C	Model tray for HSS-88	130,00
168166 C	Higher base for trimmers	150,00
320006 C	1 m water supply hose ½", with cloth inlet and 2 hose clips	25,50
320017 C	1 m drainage hose 1", translucent, made of plastic	32,50
629002 C	Sealing lubricant for trimmer, 25 ml	9,00
164026 C	Puller for grinding disc exchange for HSS-99/ HSS-88/ HSS-AZ/ HSS-ZA	46,50

Spare grinding discs

164063 C	Spare diamond grinding disc for HSS-99, entirely diamonded, directly coated, Ø 300 mm, incl. puller and sealing lubricant	2.375,00
168063 C	Spare diamond grinding disc for HSS-88, entirely diamonded, directly coated, Ø 260 mm, incl. puller and sealing lubricant	2.020,00
162063 C	Spare Dia-Quick grinding disc for HSS-AZ/-88/-ZA, partly diamonded, directly coated, Ø 255 mm, incl. sealing lubricant	1.170,00
168024 C	Spare corundum grinding disc for HSS-88, 10 mm thick, Ø 260 mm, incl. puller and sealing lubricant	363,00
160020 C	Spare corundum grinding disc for HSS-65/HSS-68, 10 mm thick, Ø 242 mm, incl. puller and sealing lubricant	398,00
162013 C	Spare foil grinding disc for HSS-AZ/-88/-ZA, 1 set = 4 pieces, Ø 250 mm, incl. sealing lubricant	182,00
162023 C	Spare foil grinding disc (fine) for HSS-AZ/-88/-ZA, 1 set = 4 pieces, Ø 250 mm, incl. sealing lubricant	220,00
168071 C	Conversion kit for HSS 88 for Dia-Quick and foil grinding disc incl. puller for grinding disc change and sealing lubricant	343,00

Modelling**Laboratory technology machines****HSM-50 / W&H Perfecta 600**

1000 – 50000 min⁻¹, 200 Watt, infinitely adjustable motor,
 engine torque 7,8 N/cm, cpl. with rest for handpiece

111986 A	Motor handpiece with foot-operated starter with blow-out function operating controls as accessories available	2.505,00
111987 A	Motor handpiece with knee-operated starter with blow-out function incl. operating controls foot-operated starter with blow-out function as accessories available	2.770,00
111988 A	Motor handpiece with table-operated starter incl. operating controls and foot-operated starter with blow-out function	2.770,00

HSM-40 / W&H Perfecta 300

1000 – 40000 min⁻¹, 160 Watt, infinitely adjustable motor,
 engine torque 5,5 N/cm, cpl. with rest for handpiece

111980 A	Motor handpiece with foot-operated starter operating controls as accessories available	1.865,00
111981 A	Motor handpiece with knee-operated starter incl. operating controls foot-operated starter as accessories available	2.165,00
111982 A	Motor handpiece with table-operated starter incl. operating controls and foot-operated starter	2.165,00

Accessories for laboratory technology machines

111121 C	Motor handpiece 50 / 40 / 300 / 600 (without supply line)	1.620,00
111117 C	Operating controls HSM-50 / W&H Perfecta 600	351,00
111111 C	Operating controls HSM-40 / W&H Perfecta 300	351,00
111124 C	Holder for operating controls 50 / 40 / 300 / 600	50,50
111119 C	Foot starter for 50* / 40 / 300 / 600*	234,00

for knee-/ and table-operated starter

*HSM-50 / W&H Perfecta 600 with blow out function

Steam cleaners
Wasi-Steam Compact and Wasi-Steam 3
(fixed water connection required)

Thorough and efficient steam cleaner suitable for all cleaning jobs in the lab. Choice of wet or dry steam

174989 B	Wasi-Steam 3 Stainless Steel	3.655,00
174987 A	<i>Wasi-Steam Compact Stainless Steel with filling container</i>	4.025,00
174995 C	Wasi-Clean Water softener (strongly recommended!) Table top unit W 190 x H 405 x D 260 mm	298,00

Wasi-Steam Classic II (manual filling)

Compact steam cleaner with steam gun, manual filling
 Flexible choice of location. Operate with distilled/softened water recommended!

174996 A	#CCWH# Wasi-Steam Classic II Stainless Steel White	1.730,00
174996 A	#CCVA# Wasi-Steam Classic II Stainless Steel	1.730,00
174970 C	<i>Wasi-Steam Box #CCWH# Stainl. Steel White, remaining stock!</i> <i>Compact working cabin with integrated suction</i>	2.080,00
174970 C	Wasi-Steam Box #CCVA# Stainless Steel Compact working cabin with integrated suction	2.080,00

Price in €
 235,00

174633 C Quantity Control Set for Wasi-Steam 3

To regulate wet steam saturation and steam quantity.

It would require to retrofit the unit at our factory or at one of our authorized service stations.

Accessories Wasi-Steam Classic/ Classic II

655002 C Lime remover, 1000 ml 105,00

174358 C Set for wall mount 23,00

390100 C Tank Closing Cap WSC 30,50

Splash protectors

174460 C Splash protection II 250,00
 W 250 x H 265 x D 255 mm
 open at the front and top
 entirely made of stainless steel

174360 C Splash protection 250,00
 W 235 x H 385 x D 294 mm
 top front with fixed stainless steel hood
 (hood: open top, chimney effect)
 with plastic basement at the bottom

174230 C Splash protection stainless steel Wasi-Steam 3 396,00

Boil-Out and polymerisation units

170970 D Wapo-Ex 12 II	9.385,00
Boil out unit for 12 flask-halves, complete with baskets, equipped with touch screen for fast and easy operation, incl. lime remover and wax remover	
170977 C Wapo-Ex 8 Plus (combination automat)	5.220,00
Compact tabletop automat for boil out up to 8 flask halves or polymerisation up to 8 flasks in a clamp, equipped with touch screen for fast and easy operation, complete with baskets and incl. service set	
170980 C Wapo-Mat III	2.405,00
Polymerisation unit for longtime- and boil polymerisation, for up to 9 flasks in the strainer, equipped with touch screen for fast and easy operation, incl. lime remover	
170993 D Polyquick 260-T105 (combination automat)	4.555,00
Pressure polymerisation unit for hot and cold polymerisation, pressure polymerisation up to 105 °C in a water bath, can also be used as a pressure compactor, very large capacity Ø 240 mm, usable depth 160 mm, suitable for extra large articulators, touch screen for fast and easy operation	
170976 C Polyquick 160-T95 (combination automat)	3.655,00
Pressure polymerisation unit for hot and cold polymerisation, pressure polymerisation up to 95 °C in a water bath can also be used as a pressure compactor, capacity Ø 150 mm, depth 285 mm touch screen for fast and easy operation Optional feature: Dry pressure compression up to approx. 70°C	

Price in €**Accessories for boil out and polymerisation units**

170120 C	Moveable stand for Wapo-Ex 8 Plus/ Wapo-Ex 8 III/ Wapo-Ex 8 with joint rollers, W x H x D: 620 x 550 x 430 mm	698,00
170100 C	Moveable stand for Wapo-Ex with joint rollers, W x H x D: 560 x 550 x 410 mm	698,00
170330 C	Service set for Wapo-Ex 12 II/ Wapo-Ex 12 consisting of lime remover and wax remover	222,00
170021 C	Service set for Wapo-Ex 8 Plus/ Wapo-Ex 8 III/ Wapo-Ex 8 consisting of 4 filter mats, lime remover, wax remover	272,00
170020 C	Service set for Wapo-Ex consisting of 10 filter mats, lime remover, wax remover	272,00
611007 C	Spare filter for Wapo-Ex 8 Plus/ Wapo-Ex 8 III/ Wapo-Ex 8 / Wapo-Ex old up to #40189 consisting of 4 filter mats 400 x 325 mm	65,50
611013 C	Spare filter for Wapo-Ex consisting of 10 filter mats 400 x 145 mm	65,50
655001 C	Wax remover , 2 liter	130,00
655002 C	Lime remover , 1000 ml	105,00

! **When you order new filters,**
please mention the S/N serial number of the unit.

!**Compactors**

170993 D	Polyquick 260-T105	4.555,00
170976 C	Polyquick 160-T95	3.655,00
180994 C	Compactor DEG-1 old item no 180999	2.495,00
380032 C	O-Ring for DEG / Polyquick 160	31,95
380063 C	O-Ring for Polyquick 260	84,00

Mixing and dosing unit

172970 C	Sidomix III	1.990,00
	Dosage and mixing unit for addition-vulcanizing duplicating silicone, 1:1	

Accessories for mixing and dosing unit

172640 C	Set material cask, cpl. 2 x 2 liter	119,00
	incl. quick fitting coupling	
172633 C	Connecting tubes cpl., 2 pcs.	77,00
172634 C	Adapter set for WaSil silicone 2 x 1 kg cask	69,50
172636 C	Adapter set for WaSil silicone 2 x 6 kg cask	77,00
622002 C	Mixing canulas 1:1, 5.4 mm inner Ø, pink, 100 pcs.	103,00
	connectable to dosage automat 1:1 Sidomix II, Sidomix, DA 2000 and similar equipment	

**! For a long durability of the mixing and dosing unit
we recommend the use of our Wassermann
WaSil duplicating silicone. More on Page 26.**

!

Ejection press

128999 A APR-1	1.180,00
for cast rings/ cast flasks (DeguDent/ Heraeus Kulzer), cast ring set 1	

Flask presses

126999 A Hydraulic press WW-33 for 1-3 flasks	1.520,00
126998 A Hydraulic press WW-33 for 2-4 flasks	1.640,00
129997 C Universal flask clamp for 1 flask	221,00
129996 C Universal flask clamp for 2 flasks	229,00
129995 C Universal flask clamp for 3 flasks	234,00

Accessories for flask presses

129994 C Stainless steel clamp clip for 1 flask	166,00
129993 C Stainless steel clamp clip for 2 flasks	174,00
129992 C Stainless steel clamp clip for 3 flasks	184,00
129990 C Universal flask	144,00
629005 C Spare bottle hydraulic oil, 500 ml	105,00
for WW-11, WW-22, WW-33 and KÜVOMAT	
629004 C Hydraulic Oil, Bottle 100 ml	37,50

Surface finish**Fine blasting units**

190986 A	Cemat-4 Plus with intuitive hand controls with 4 blasting material containers and sand-drawer, individually adjustable working pressure per chamber	7.655,00
190984 A	Cemat-4 with 4 blasting material containers and sand-drawer, individually adjustable working pressure per chamber	6.570,00
190985 A	Cemat-3 Plus with intuitive hand controls with 3 blasting material containers and sand-drawer, individually adjustable working pressure per chamber	6.430,00
190983 A	Cemat-3 with 3 blasting material containers and sand-drawer, individually adjustable working pressure per chamber	5.900,00
190981 A	Cemat-2 Plus with intuitive hand controls with 2 blasting material containers and sand-drawer, individually adjustable working pressure per chamber	5.780,00
190982 A	Cemat-2 with 2 blasting material containers and sand-drawer, individually adjustable working pressure per chamber	5.250,00
190884 A	Gobi Comfort-4 with 4 blasting material containers and sand-drawer,	4.910,00
190883 A	Gobi Comfort-3 with 3 blasting material containers and sand-drawer,	4.285,00
190882 A	Gobi Comfort-2 with 2 blasting material containers and sand-drawer,	3.605,00
190881 A	Gobi Classic with 2 blasting material containers	2.505,00
190885 A	Puna II with 2 blasting material containers	1.730,00

! The Wassermann fine blasting units are optimally designed for use with our Wassermann dental extraction systems. Page 22.

!

For optimal durability, we recommend Wassermann blasting material to extend the service life of our units. Page 27.

Accessories for fine blasting units**Accessories Cemat II/NT, Gobi Comfort/Classic, Puna/II**

For Cemat OLD please see page 20

190553 C	Spare fine blasting nozzle 0.4 mm (only Cemat-NT)	112,00
190550 C	Spare fine blasting nozzle 0.6 mm	112,00
190551 C	Spare fine blasting nozzle 0.8 mm	112,00
190552 C	Spare fine blasting nozzle 1.2 mm	112,00
190565 C	Spare nozzle (tank) 2.0 mm	19,00
190566 C	Spare nozzle (tank) 1.4 mm	19,00
190494 C	Key for assembly/ disassembly of the tank nozzle	34,50
190533 C	Rep. Set Sleeve Gobi Cemat II	74,00

Accessories for Cemat II

190489 C	Spare hose -yellow-	111,00
190485 C	Spare hose -blue-	111,00
190475 C	Spare hose -red-	111,00
190487 C	Spare hose -green-	111,00
190432 C	Spare protection sheet with special coating	112,00

Accessories for Gobi Comfort and Gobi Classic

190739 C	Spare hose -yellow-	105,00
190738 C	Spare hose -blue-	105,00
190729 C	Spare hose -red-	105,00
190728 C	Spare hose -green-	105,00
190432 C	Spare protection sheet with special coating	112,00

Accessories for Puna II

190475 C	Spare hose -red-	111,00
190485 C	Spare hose -blue-	111,00
190686 C	Spare protection sheet with special coating	112,00

Accessories for fine blasting units (previous versions)**Accessories for Cemat-NT****(Item no. 190992/ 190990/ 190993)**

Spare Nozzles see page 18

190432 C	Spare protection sheet with special coating from serial number 128538	112,00
190458 C	Conversion kit protection sheet with special coating up to serial number 128537	129,00
190230 C	Spare protection sheet up to serial number 128537	112,00
612003 C	Spare protection foil (plastic), 1 set = 5 pieces, up to serial number 128537	85,50
190239 C	Spare hoses	129,00
190533 C	Rep. Set Sleeve Gobi Cemat II	74,00

Accessories for Gobi-2 (Item no. 190896)

Spare Nozzles see page 18

190432 C	Spare protection sheet with special coating from serial number 128423	112,00
190636 C	Conversion-kit protection sheet with special coating # 124048 - # 127319	405,00
612004 C	Spare protection foil (plastic), 1 set = 5 pieces, up to serial number 127319	85,50
190487 C	Spare hose for green chamber	111,00
190475 C	Spare hose for red chamber	111,00
190485 C	Spare hose for blue chamber	111,00
190489 C	Spare hose for yellow chamber	111,00
190533 C	Rep. Set Sleeve Gobi Cemat II	74,00

Accessories for Puna (Item no. 190897)

Spare Nozzles see page 18

190686 C	Spare protection sheet with special coating	112,00
190475 C	Spare hose for red chamber	111,00
190485 C	Spare hose for blue chamber	111,00
190533 C	Rep. Set Sleeve Gobi Cemat II	74,00

Polishing compact units

146979 B	Polishing unit WP-Ex 3000 II incl. adjustable motor (750-6000 min ⁻¹), dust extraction, LED illumination and LED Spot, rubber polishing troughs, drawer, brush spindle left + right, filters, 2 storage bowls	5.845,00
146978 B	Polishing unit WP-Ex 2000 II incl. 500 W motor, dust extraction, LED illumination, rubber polishing troughs, drawer, brush spindle left + right, filters, 2 storage bowls	4.760,00
146971 B	Polishing unit WP-Ex 10 II incl. 500 W motor, dust extraction, LED illumination, rubber polishing trough, brush spindle right , filter, 1 storage bowl	3.125,00

Polishing motor

155999 B	WSM-2 (1500/ 3000 min⁻¹) incl. brush spindle left + right	1.125,00
-----------------	--	----------

Polishing troughs

143123 C	Polishing trough without extractor nozzle incl. adjustable protection sheet and removable rubber trough	271,00
143115 C	Polishing trough with extractor nozzle (Ø75 mm), incl. cap incl. adjustable protection sheet and removable rubber trough	305,00
143131 C	Polishing trough with extractor nozzle (Ø40 mm), incl. cap incl. adjustable protection sheet and removable rubber trough	305,00
143550 C	Extraction traverse II, with switch to guide the air output nozzle Ø75 mm	608,00

Accessories for polishing lathes

152025 C	Brush spindle 92 mm for left side	56,50
152022 C	Brush spindle 92 mm for right side	56,50
152057 C	Stone holder for left side	56,50
152056 C	Stone holder for right side	56,50
152044 C	Brush spindle 125 mm for left side	85,50
152043 C	Brush spindle 125 mm for right side	85,50
152049 C	Chuck clamping from 1-6 mm	222,00
155123 C	Wrench for pull nut	69,00
146390 C	Daylight LED Spot	307,00
146585 C	Conversion Kit with two glass panels and axle from old version to current version	479,00
146588 C	RepSet Glass Panel	224,00

Filter for polishing compact units

611005 C	Spare filter bag for WP-Ex 3000 II/ 2000 II/ 10 II dust class M, for micro fine dust, 1 set = 2 pieces	100,00
611011 C	Spare filter bag WP-Ex 3000/ 2000/ 10 dust class M, for micro fine dust, 1 set = 2 pieces	85,50
611021 C	Spare filter bag for WP-Ex 80 dust class M, for micro fine dust, 1 set = 2 pieces	85,50

! Please sent us the serial number of your unit in case of ordering filter material. The S/N label is located on the back or right side of your unit. !

Automatic stand

146170 D	Automatic stand for WP-Ex 3000 II/ 2000 II vertically adjustable from 74 – 109 cm Dimension (table top): D 560 x W 675 mm	3.125,00
-----------------	--	----------

Price in €**Suction technology**

140997 A	Single-place-dust extraction SG-1/2 D (Suction port 75 mm) Brushless motor. Equipped with automatic switch-off when the filter is full. Filter included. Without accessories.	1.920,00
140995 A	Single place dust extraction SG-1/1 (Suction port 75 mm) Brushless motor. Filter included. Without accessories.	1.155,00
140991 C	Single place dust extraction SG-1/1 Comfort (Suction port 75 mm) Brushless motor. With automatic fill level monitoring. Equipped with auto-start-system.	1.455,00
140970 C	Suction unit SG-Economy (Suction port 40 mm) incl. 3 stage filter system, without accessories. Recommended for HSS-TR (dry trimmer) Powerful brushed motor for reliable performance. Equipped with auto-start-system.	1.410,00

Price in €

145997 A	Dust extraction drawer LSG-02 DE complete with universal mounting frame and automatic switch-off in case of exhausted filter capacity, incl. filter , additional equipped with auto-start-system	2.435,00
145998 A	Dust extraction drawer LSG-02 D complete with universal mounting frame and automatic switch-off in case of exhausted filter capacity, incl. filter	2.245,00
145240 C	Activated charcoal filters for monomer for LSG-02 DE/ LSG-02 D for elimination of dangerous vapors	715,00
145200 C	Universal extractor hood for LSG-02 D and DE with exchangeable protection sheet	428,00
145005 C	Additional drawer for LSG-02 D and DE can also be retrofitted	232,00
149999 A	Mobile dust extraction Timobil cpl. with handpiece motor 1.000 – 40.000 min ⁻¹ , 160 Watt, infinitely adjustable, standard height 900 mm	4.760,00
149140 C	Extractor hood illuminated for Timobil	641,00

Suction Box

143372 C	Airbox One New! Compact. Mobile. Powerful. *Mobile extraction box and extraction unit in one device *With practical carrying handle for flexible use *Incl. daylight LED lighting *Powerful extraction performance *M-class filter	1.390,00
143390 C	Suction Box Compact III Incl. daylight LED lighting With 1 external extraction port Ø 75 mm Extraction port can be mounted at the front or rear	820,00

Accessories for Suction Box Compact III +previous versions

143303 C	Set of arm rests	87,50
143305 C	Set of cover plates	225,00
190686 C	Spare protection sheet with special coating	112,00

Price in €

Accessories for suction units**Telescopic channels + corresponding suction parts**

143070 C	Telescopic channel 570 mm	155,00
143090 C	Telescopic channel 455 mm	165,00
143074 C	Telescopic channel 480 mm	164,00
143066 C	Suction part, transverse	123,00
143080 C	Suction part, straight	135,00
145245 C	Activated charcoal filters for monomer - suction part, straight for elimination of dangerous vapors	715,00
143108 C	Cap for suction part transverse	51,00
143200 C	Universal hood for suction part straight	428,00
143220 C	Suction part with wooden rest	357,00
143221 C	Suction part with wooden rest/plastic screen	430,00
143231 C	Suction part with wooden rest/glass screen	479,00

Special channel + appropriate suction parts

(e.g. for installation in Freuding tables)

143093 C	Special channel with hose coupling	187,00
143041 C	Extractor nozzle	225,00

! Channel and one suction part are always both required!

!**Hoses**

140057 C	PVC-hose 75 Ø - 1.20 m length, incl. 2 hose clips	112,00
320011 C	PVC-hose 75 Ø - 1.80 m length	159,00
320010 C	PVC-hose 75 Ø - 2.40 m length	196,00
320012 C	PVC-hose 75 Ø - 3.60 m length	274,00
322003 C	Hose clip 75 Ø	5,40
327023 C	PVC-hose 51 Ø inner - per meter	61,50

Special accessories

140056 C	Throttle valve 75 Ø	171,00
-----------------	---------------------	--------

! Reduction- and transition pieces of all measurements available.

!

Would you please be so kind as to mention inner- and outer diameter in case of supplementary order.

Price in €**Filters for suction units****Filters & accessories for suction unit SG-10**

611010 C	Spare filter bag yellow dust class M, for micro fine dust, 1 set = 4 pieces	200,00
611018 C	Spare filter bag white, for plaster class 1 and 2 only (!) regenerative, 1 set = 4 pieces	200,00
610073 C	*Filter bag white (plaster class 1/ class 2), big	79,00
610074 C	*Filter bag yellow (micro fine dust), big	79,00
<i>* to be used only in combination with item no. 140410</i>		
140410 C	Filter extension module SG-10	276,00
140435 C	Filter cassette, cpl. SG-10/ SG-20	93,00
140424 C	Module activated carbon filter SG-10	549,00
140436 C	Activated carbon filter (AKM 6-ply) SG-10	157,00

SG-1/1 + SG-2/2 + Timobil

611010 C	Spare filter bag yellow, dust class M, for micro fine dust, 1 set = 4 pieces	200,00
611018 C	Spare filter bag white, (class 1/class 2) for plaster, regenerable, 1 set = 4 pieces	200,00

SG-1/2 D+ LSG-02 DE/ LSG-02 D

611020 C	Spare filter bag yellow, dust class M, for micro fine dust, 1 set = 4 pieces	240,00
611017 C	Spare filter bag white, (class 1/class 2) for plaster, regenerable, 1 set = 3 pieces	182,00
610056 C	Spare filter monomer for LSG-02	72,50
611022 C	<u>Filter-Kit SG-Economy</u> 1 x Spare-main filter / cartridge filter SG-Economy 5 x Spare-filter bag SG-Economy 2 x Spare-HEPA exhaust filter SG-Economy	236,00

Filters & accessories for suction unit SG-20

610073 C	Filter bag white (plaster class 1/ class 2)	79,00
610074 C	Filter bag yellow (micro fine dust)	79,00
140551 C	Activated carbon filter, cpl. SG-20	152,00
140435 C	Filter cassette, cpl. SG-10/ SG-20	93,00
610019 C	Motor protective filter SG-20	12,50

WTS-1R

611019 C	filter bag FB 25, 1 Set = 3 pieces	112,00
-----------------	------------------------------------	--------

Price in €

Consumable materials**Duplicating silicone****WaSil Soft**

Addition-vulcanizing duplicating silicone, 1:1
hardness Shore A: 16-18
component A white (catalyst)
component B **yellow** (base)

656060 C	WaSil soft 2 x 1 kg component A/ B	119,00
656061 C	WaSil soft 2 x 6 kg component A/ B	605,00
656062 C	WaSil soft 2 x 25 kg component A/ B	2.395,00

WaSil Speed

Addition-vulcanizing duplicating silicone, 1:1
high density, hardness Shore A: > 24
component A white (catalyst)
component B **orange** (base)

656080 C	WASIL SPEED A/B 1 Kg	135,00
656081 C	WASIL SPEED A/B 6 Kg	745,00
656082 C	WASIL SPEED A/B 25 Kg	2.835,00

Mixing canulas, 1:1
5.4 mm inner Ø

622002 C	mixing canulas, pink, 100 pcs.	103,00
	connectable to dosage automat 1:1 Sidomix II, Sidomix, DA 2000 and similar equipment	

Price in €**Blasting materials****Aluminium oxyde**

For blasting of all dental alloys. The purity, sharp-edged grains and extreme hardness guarantee very efficient blasting.

Silicosis-free, made in Germany

656010 C	Aluminium oxyde 50 µm – 10 kg	97,50
656015 C	Aluminium oxyde 50 µm – 20 kg	201,00
656013 C	Aluminium oxyde 110 µm - 10 kg	97,50
656016 C	Aluminium oxyde 110 µm - 20 kg	201,00
656014 C	Aluminium oxyde 250 µm - 10 kg	97,50
656018 C	Aluminium oxyde 250 µm – 20 kg	201,00

Glass beads

For gentle cleaning and compacting of sensitive surfaces of precious and non-precious alloys. For shining and smoothing out of crowns, bridges and CrCo techniques.

Silicosis-free, made in Germany

656020 C	Glass beads 50 µm – 8 kg	76,00
656021 C	Glass beads 50 µm – 15 kg	107,00
656019 C	Glass beads 125 µm – 8 kg	76,00
656022 C	Glass beads 125 µm – 15 kg	107,00



Wassermann Dental-Maschinen GmbH
Rudorffweg 15–17 · 21031 Hamburg · Germany
Phone: +49 (0)40 730 926-0
export@wassermann-dental.com · www.wassermann.hamburg



[/WassermannDentalmaschinen](#)