

Recommended Retail Prices 2024



## **PRICE LIST 2024**

This price list is effective from 01.01.2024, ex works Hamburg and is subject to change without prior notice and cannot be considered binding to the seller for any specific period.

Prices are quoted in EUR (EURO), without V.A.T. and subject to our normal terms and conditions of sale regarding delivery and payment. You could request these terms by phone or download them from our homepage.

Specifications are subject to alteration without prior notice.

Orders placed below EUR 50.-- are subject to a handling fee of EUR 10.--.

The plastic coating usually is:  
natural white close to RAL 9016 and black close to RAL 9005.

Colours are also subject to change without notice. Alternative colours are available at additional surcharge.

All units are available in 220-240 V / 50/60 Hz.  
Other voltages upon enquiry. Please let us know.  
Units equipped with heating are only available in 220-240 V / 50/60 Hz.

**Wassermann Dental-Maschinen GmbH**  
**Rudorffweg 15 – 17, 21031 Hamburg, Germany**

Export Department:

Phone: +49 40 730 926 22 or +49 40 730 926 24

[export@wassermann-dental.com](mailto:export@wassermann-dental.com)

[www.wassermann.hamburg](http://www.wassermann.hamburg)

[www.facebook.com/WassermannDentalmaschinen](https://www.facebook.com/WassermannDentalmaschinen)

## Contents

Dental Cleaner.....	4
Modell production.....	4
<a href="#">Plaster Dispenser and Accessories.....</a>	<a href="#">4</a>
<a href="#">Dosing scale.....</a>	<a href="#">4</a>
<a href="#">Accessories for DSW-2/ -2/4K.....</a>	<a href="#">4</a>
<a href="#">Vacuum mixers.....</a>	<a href="#">5</a>
<a href="#">Accessories for vacuum mixers.....</a>	<a href="#">5</a>
<a href="#">Vibrators.....</a>	<a href="#">6</a>
<a href="#">Accessories for vibrators.....</a>	<a href="#">6</a>
<a href="#">Trimmers.....</a>	<a href="#">7</a>
<a href="#">Accessories for trimmers.....</a>	<a href="#">8</a>
<a href="#">Spare grinding discs.....</a>	<a href="#">9</a>
Modelling.....	10
<a href="#">Laboratory technology machines.....</a>	<a href="#">10</a>
<a href="#">Accessories for laboratory technology machines.....</a>	<a href="#">10</a>
<a href="#">Steam cleaners.....</a>	<a href="#">11</a>
<a href="#">Accessories Wasi-Steam Compact.....</a>	<a href="#">11</a>
<a href="#">Accessories Wasi-Steam Classic/ Classic II.....</a>	<a href="#">12</a>
<a href="#">Splash protection.....</a>	<a href="#">12</a>
<a href="#">Boil out and polymerisation units.....</a>	<a href="#">13</a>
<a href="#">Accessories for boil out and polymerisation units.....</a>	<a href="#">14</a>
<a href="#">Compactors.....</a>	<a href="#">14</a>
<a href="#">Mixing and dosing unit.....</a>	<a href="#">15</a>
<a href="#">Accessories for mixing and dosing unit.....</a>	<a href="#">15</a>
<a href="#">Ejection press.....</a>	<a href="#">16</a>
<a href="#">Flask presses.....</a>	<a href="#">16</a>
<a href="#">Accessories for flask presses.....</a>	<a href="#">16</a>
Surface finish.....	17
<a href="#">Fine blasting units.....</a>	<a href="#">17</a>
<a href="#">Accessories Cemat II, Gobi Comfort, Gobi Classic, Puna II.....</a>	<a href="#">18</a>
<a href="#">Accessories for Cemat II.....</a>	<a href="#">18</a>
<a href="#">Accessories for Gobi Comfort and Gobi Classic.....</a>	<a href="#">18</a>
<a href="#">Accessories for Puna II.....</a>	<a href="#">18</a>
<a href="#">Accessories for Cemat-NT (Item no. 190992/ 190990/ 190993).....</a>	<a href="#">19</a>
<a href="#">Accessories for Gobi-2 (Item no. 190896).....</a>	<a href="#">19</a>
<a href="#">Accessories for Puna (Item no. 190897).....</a>	<a href="#">19</a>
<a href="#">Accessories for Cemat (Item no. 190996/ 190999/ 190998/ 190997).....</a>	<a href="#">20</a>
<a href="#">Polishing compact units.....</a>	<a href="#">20</a>
<a href="#">Polishing motor.....</a>	<a href="#">20</a>
<a href="#">Polishing troughs.....</a>	<a href="#">20</a>
<a href="#">Accessories for polishing compact units/ polishing motors.....</a>	<a href="#">21</a>
<a href="#">Filter for polishing compact units.....</a>	<a href="#">21</a>
<a href="#">Automatic stand.....</a>	<a href="#">21</a>
Suction technology.....	22
<a href="#">Suction Box.....</a>	<a href="#">23</a>
<a href="#">Accessories for suction units.....</a>	<a href="#">24</a>
<a href="#">Filters for suction units.....</a>	<a href="#">25</a>
Consumable materials.....	26
<a href="#">Duplicating silicone.....</a>	<a href="#">26</a>
<a href="#">Mixing canulas, 1:1.....</a>	<a href="#">26</a>
<a href="#">Blasting materials.....</a>	<a href="#">27</a>

Price in €

**Modell production****Plaster Dispenser and Accessories**

<b>166997 D</b>	<b>Plaster Dispenser WGS-25W II</b> Wall Model, rotatable (180 °), for 25 kg plaster container	4.460,00
<b>166115 D</b>	<b>Table-Top Support</b> for WGS-25 W II	658,00

**Dosing scales**

<b>172988 A</b>	<b>DSW-2 (with 1 pump, not cooled)</b> For exact dosage of plaster and investment material You can order up to 5 additional pumps	2.399,00
-----------------	---	----------

**Accessories for DSW-2/ -2/4K**

<b>172299 C</b>	Additional pump equipment (not cooled)	218,00
-----------------	--	--------

Price in €

**Vacuum mixers**

<b>171970 A</b>	<b>Wamix-Touch</b> with vacuum pump, incl. wall holder and mixing cup 350 ml with touch screen for fast and easy operation	2.149,00
<b>171971 A</b>	<b>Wamix-Touch Injector</b> with injector technique, incl. wall holder and mixing cup 350 ml with touch screen for fast and easy operation	1.939,00
<b>171972 B</b>	<b>Wamix-Economy <span style="color: red;">New!</span></b> with vacuum pump	1.399,00

**Accessories for vacuum mixers****for Wamix-Touch/ Wamix-Touch Injector and  
Wamix-2/ -2M/ -2MS/ -3M/ -2000/ CCI<sup>1</sup>**CCI<sup>1</sup> = Wamix-Comfort Plus, -Comfort, -Classic, -Injector

<b>171295 C</b>	Mixing cup Wamix mini, 75 ml (75 ml mixing cup not suitable for beaker mount Wamix-Comfort Plus)	209,00
<b>171230 C</b>	Mixing cup Wamix small, 200 ml	218,00
<b>171283 C</b>	Mixing cup Wamix medium, 350 ml	225,00
<b>171237 C</b>	Mixing cup Wamix large, 600 ml	233,00
<b>171269 C</b>	Mixing cup Wamix extra large, 1000 ml (1000 ml mixing cup not suitable for beaker mount Wamix-Comfort Plus)	297,00
<b>171744 C</b>	Table stand Wamix-Touch/ Wamix-Touch Injector	341,00
<b>171777 C</b>	Table stand Wamix-CCI <sup>1</sup> for model IV and model 2012	341,00
<b>171760 C</b>	Table stand Wamix for models before 2012	319,00
<b>171122 C</b>	Set of sealings (set = 5 pcs.)	8,50
<b>610025 C</b>	Sinterfilter	23,00

<b><u>Vibrators</u></b>		<b>Price in €</b>
<b>169997 A</b>	<b>KV-56, <u>plastic-powder-coated</u> cast housing</b> White Table dimensions 430 x 286 mm, vibration intensity infinitely variable	1.609,00
<b>169999 A</b>	<b>KV-36, <u>plastic-powder-coated</u> stainless steel</b> White Table dimensions 293 x 200 mm, vibration intensity infinitely variable	759,00
<b>169992 A</b>	<b>KV-26 plus, <u>plastic-powder-coated</u> stainless steel</b> White Table dimensions 210 x 130 mm, vibration intensity infinitely variable	657,00
<b>169996 A</b>	<b>KV-26, <u>plastic-powder-coated</u> stainless steel</b> White Table dimensions 210 x 130 mm	543,00
<b>169998 A</b>	<b>KV-16, <u>plastic-powder-coated</u> stainless steel</b> White Table dimensions 130 x 130 mm	469,00
<b>169989 A</b>	<b>KV-36, <u>matt</u> stainless steel</b> Table dimensions 293 x 200 mm, vibration intensity infinitely variable	779,00
<b>169988 A</b>	<b>KV-26 plus, <u>matt</u> stainless steel</b> Table dimensions 210 x 130 mm, vibration intensity infinitely variable	666,00
<b>169987 A</b>	<b>KV-26, <u>matt</u> stainless steel</b> Table dimensions 210 x 130 mm	550,00
<b>169986 A</b>	<b>KV-16, <u>matt</u> stainless steel</b> Table dimensions 130 x 130 mm	470,00
 <b><u>Accessories for vibrators</u></b>		
<b>169510 C</b>	Spare rubber plate KV-56	132,00
<b>169010 C</b>	Spare rubber plate KV-36	62,00
<b>169101 C</b>	Spare rubber plate KV-26/ KV-26 plus	54,00
<b>169201 C</b>	Spare rubber plate KV-16	30,50

	<u>Trimmers</u>	Price in €
164998 A	<b>HSS-99 - 1200 Watt – incl. diamond grinding disc</b> grinding disc entirely diamonded, directly coated, incl. solenoid valve and ball tap, large grinding area	4.765,00
168998 A	<b>HSS-88 - 1000 Watt - incl. diamond grinding disc</b> incl. solenoid valve, grinding disc entirely diamonded, directly coated	3.750,00
168996 A	<b>HSS-88 - 1000 Watt - incl. Dia-Quick grinding disc</b> incl. equipment for first installation and solenoid valve, grinding disc partly diamonded, directly coated	3.299,00
163990 B	<b>HSS-TR Dry Trimmer Dia-Quick <span style="color: red;">New!</span></b> with partially diamond disc	3.155,00
168999 A	<b>HSS-88 - 1000 Watt – incl. corundum-grinding disc</b> incl. solenoid valve, corundum-grinding disc 10 mm thick	2.285,00
168991 A	<b>HSS-88 - 1000 Watt - incl. foil grinding disc</b> incl. solenoid valve, incl. equipment for first installation	2.285,00
168997 A	<b>HSS-88 - 1000 Watt – without grinding disc</b> incl. solenoid valve <b>all WASSERMANN-grinding discs are applicable</b>	2.040,00
162997 A	<b>HSS-AZ - 750 Watt – incl. Dia-Quick grinding disc</b> incl. solenoid valve, grinding disc partly diamonded, directly coated	2.549,00
162996 A	<b>HSS-AZ - 750 Watt - incl. foil grinding disc</b> incl. solenoid valve	1.599,00
161996 A	<b>HSS-ZA</b> <b>incl. Dia-Quick grinding disc and brush spindle on left side</b> combination trimmer/polishing motor -500 Watt-	3.575,00
161997 A	<b>HSS-ZA</b> <b>incl. foil grinding disc and brush spindle on left side</b> combination trimmer/polishing motor -500 Watt-	2.490,00
161114 C	<b>Polishing trough base for HSS-ZA (2017)</b> (Polishing troughs see page 19)	76,50

<b><u>Accessories for trimmers</u></b>		<b>Price in €</b>
<b>Grinding attachment for HSS-99 for producing 3D orthodontic models:</b>		
<b>164130 D</b>	Model table KMS according to Begg	2.175,00
<b>164110 D</b>	Laser light guide for model table	828,00
<b>164034 D</b>	Orthodontic template, lower jaw/ upper jaw/ tuber plane	51,00
<b>Grinding attachment for HSS-88 for producing 3D orthodontic models:</b>		
<b>168140 D</b>	Model table KMS according to Begg	2.175,00
<b>168120 D</b>	Laser light guide for model table	911,00
<b>164034 D</b>	Orthodontic template, lower jaw/ upper jaw/ tuber plane	51,00
<b>!</b>	<b>The Grinding attachment KMS can be delivered also as Swiss and EOS-versions. Price on demand.</b>	<b>!</b>
<b>164161 C</b>	Tray for trimmer HSS-99 fitting in a sink (We always need dimensions of the sink and design 1:1)	313,00
<b>168160 C</b>	Tray for trimmer HSS-88/ AZ fitting in a sink (We always need dimension of the sink and design 1:1)	270,00
<b>164037 C</b>	Inclined bearing table for HSS-99	102,00
<b>168037 C</b>	Inclined bearing table for HSS-88/HSS-AZ/HSS-ZA	100,00
<b>168168 C</b>	Model tray for HSS-88	123,00
<b>168166 C</b>	Higher base for trimmers	140,00
<b>320006 C</b>	1 m water supply hose ½", with cloth inlet and 2 hose clips	24,00
<b>320017 C</b>	1 m drainage hose 1", translucent, made of plastic	30,50
<b>710001 C</b>	Dressing tool for roughening of corundum grinding discs	98,50
<b>0</b>		
<b>629002 C</b>	Sealing lubricant for trimmer, 25 ml	8,50
<b>164026 C</b>	Puller for grinding disc exchange for HSS-99/ HSS-88/ HSS-AZ/ HSS-ZA	44,00



<u>Spare grinding discs</u>		Price in €
<b>164063 C</b>	<b>Spare diamond grinding disc for HSS-99,</b> entirely diamonded, directly coated, Ø 300 mm, incl. puller and sealing lubricant	2.270,00
<b>168063 C</b>	<b>Spare diamond grinding disc for HSS-88,</b> entirely diamonded, directly coated, Ø 260 mm, incl. puller and sealing lubricant	1.930,00
<b>162063 C</b>	<b>Spare Dia-Quick grinding disc for HSS-AZ/-88/-ZA,</b> partly diamonded, directly coated, Ø 255 mm, incl. sealing lubricant	1.115,00
<b>168024 C</b>	<b>Spare corundum grinding disc for HSS-88,</b> 10 mm thick, Ø 260 mm, incl. puller and sealing lubricant	347,00
<b>160020 C</b>	<b>Spare corundum grinding disc for HSS-65/HSS-68,</b> 10 mm thick, Ø 242 mm, incl. puller and sealing lubricant	380,00
<b>162013 C</b>	<b>Spare foil grinding disc for HSS-AZ/-88/-ZA,</b> 1 set = 4 pieces, Ø 250 mm, incl. sealing lubricant	173,00
<b>162023 C</b>	<b>Spare foil grinding disc (fine) for HSS-AZ/-88/-ZA,</b> 1 set = 4 pieces, Ø 250 mm, incl. sealing lubricant	210,00
<b>168071 C</b>	<b>Conversion kit for HSS 88</b> for Dia-Quick and foil grinding disc incl. puller for grinding disc change and sealing lubricant	327,00

**Modelling****Price in €****Laboratory technology machines****HSM-50 / W&H Perfecta 600**

1000 – 50000 min<sup>-1</sup>, 200 Watt, infinitely adjustable motor,  
 engine torque 7,8 N/cm, cpl. with rest for handpiece

<b>111986 A</b>	<b>Motor handpiece with foot-operated starter</b> with blow-out function operating controls as accessories available	2.399,00
<b>111987 A</b>	<b>Motor handpiece with knee-operated starter</b> with blow-out function incl. operating controls foot-operated starter with blow-out function as accessories available	2.650,00
<b>111988 A</b>	<b>Motor handpiece with table-operated starter</b> incl. operating controls and foot-operated starter with blow-out function	2.650,00

**HSM-40 / W&H Perfecta 300**

1000 – 40000 min<sup>-1</sup>, 160 Watt, infinitely adjustable motor,  
 engine torque 5,5 N/cm, cpl. with rest for handpiece

<b>111980 A</b>	<b>Motor handpiece with foot-operated starter</b> operating controls as accessories available	1.785,00
<b>111981 A</b>	<b>Motor handpiece with knee-operated starter</b> incl. operating controls foot-operated starter as accessories available	2.069,00
<b>111982 A</b>	<b>Motor handpiece with table-operated starter</b> incl. operating controls and foot-operated starter	2.069,00

**Accessories for laboratory technology machines**

<b>111121 C</b>	<b>Motor handpiece 50 / 40 / 300 / 600 (without supply line)</b>	1.550,00
<b>111117 C</b>	<b>Operating controls HSM-50 / W&amp;H Perfecta 600</b>	335,00
<b>111111 C</b>	<b>Operating controls HSM-40 / W&amp;H Perfecta 300</b>	335,00
<b>111124 C</b>	<b>Holder for operating controls 50 / 40 / 300 / 600</b>	48,00
<b>111119 C</b>	<b>Foot starter for 50* / 40 / 300 / 600*</b> for knee-/ and table-operated starter *HSM-50 / W&H Perfecta 600 with blow out function	223,00

<u>Steam cleaners</u>		Price in €
<b>Wasi-Steam Compact (fixed water connection required)</b>		
Thorough and efficient steam cleaner		
Suitable for all cleaning jobs in the lab		
Choice of wet or dry steam		
Basic version: Tabletop unit, <u>permanent water supply</u>		
<b>174994 A</b>	<b>#CCWH# Wasi-Steam Compact Stainless Steel White</b>	3.499,00
<b>174994 A</b>	<b>#CCVA# Wasi-Steam Compact Stainless Steel</b>	3.499,00
<b>174995 C</b>	<b>Wasi-Clean Water softener (strongly recommended!)</b>	278,00
	Table top unit W 190 x H 405 x D 260 mm	
<b>Wasi-Steam Classic II (manual filling)</b>		
Compact steam cleaner with steam gun, <u>manual filling</u>		
Flexible choice of location		
Operation with distilled / softened water is possible		
Fastening set for wall mounting available as accessory		
<b>174996 A</b>	<b>#CCWH# Wasi-Steam Classic II Stainless Steel White</b>	1.650,00
<b>174996 A</b>	<b>#CCVA# Wasi-Steam Classic II Stainless Steel</b>	1.650,00
<b>174970 C</b>	<b>Wasi-Steam Box #CCWH# <span style="color: red;">New!</span></b>	1.990,00
	Compact working cabin with integrated suction	

	<b>Price in €</b>
<b>174414 C Quantity Control Accessory Set</b> Wasi-Steam Compact Accessory to regulate <u>wet steam saturation</u> and <u>steam quantity</u> (optional). Retrofitting is possible, the unit can be upgraded at any time. It would require to retrofit the unit at our factory or at one of our authorized service stations.	223,00
<b>174430 C Filling container, 5 Liter, for Wasi-Steam Compact</b> Wasi-Steam Compact accessory for <u>manual filling</u> (optional) Flexible choice of location Operation with distilled / softened water is possible	367,00

#### Accessories Wasi-Steam Classic/ Classic II

<b>655002 C Lime remover, 1000 ml</b>	100,00
<b>174358 C Set for wall mount</b>	21,50
<b>390100 C Tank Closing Cap WSC</b>	28,50

#### Splash protectors

<b>174460 C Splash protection II</b> W 250 x H 265 x D 255 mm open at the front and top entirely made of stainless steel	233,00
<b>174360 C Splash protection</b> W 235 x H 385 x D 294 mm top front with fixed stainless steel hood (hood: open top, chimney effect) with plastic basement at the bottom	233,00

<b><u>Boil-Out and polymerisation units</u></b>		<b>Price in €</b>
<b>170970 D</b>	<b>Wapo-Ex 12 II</b> Boil out unit for 12 flask-halves, complete with baskets, equipped with touch screen for fast and easy operation, incl. lime remover and wax remover	8.990,00
<b>170977 C</b>	<b>Wapo-Ex 8 Plus</b> (combination automat) Compact tabletop automat for boil out up to 8 flask halves or polymerisation up to 8 flasks in a clamp, equipped with touch screen for fast and easy operation, complete with baskets and incl. service set	4.999,00
<b>170980 C</b>	<b>Wapo-Mat III</b> Polymerisation unit for longtime- and boil polymerisation, for up to 9 flasks in the strainer, equipped with touch screen for fast and easy operation, incl. lime remover	2.299,00
<b>170993 C</b>	<b>Polyquick 260-T105</b> (combination automat) Pressure polymerisation unit for hot and cold polymerisation, pressure polymerisation up to 105 °C in a water bath, can also be used as a pressure compactor, very large capacity Ø 240 mm, usable depth 160 mm, suitable for extra large articulators, touch screen for fast and easy operation	4.359,00
<b>170976 C</b>	<b>Polyquick 160-T95</b> (combination automat) Pressure polymerisation unit for hot and cold polymerisation, pressure polymerisation up to 95 °C in a water bath can also be used as a pressure compactor, capacity Ø 150 mm, depth 285 mm touch screen for fast and easy operation Optional feature: Dry pressure compression up to approx. 70°C	3.499,00

<b><u>Accessories for boil out and polymerisation units</u></b>		<b>Price in €</b>
<b>170120 C</b>	<b>Moveable stand</b> for Wapo-Ex 8 Plus/ Wapo-Ex 8 III/ Wapo-Ex 8 with joint rollers, W x H x D: 620 x 550 x 430 mm	652,00
<b>170100 C</b>	<b>Moveable stand</b> for Wapo-Ex with joint rollers, W x H x D: 560 x 550 x 410 mm	652,00
<b>170330 C</b>	<b>Service set for Wapo-Ex 12 II/ Wapo-Ex 12</b> consisting of lime remover and wax remover	212,00
<b>170021 C</b>	<b>Service set for Wapo-Ex 8 Plus/ Wapo-Ex 8 III/ Wapo-Ex 8</b> consisting of 4 filter mats, lime remover, wax remover	260,00
<b>170020 C</b>	<b>Service set for Wapo-Ex</b> consisting of 10 filter mats, lime remover, wax remover	260,00
<b>611007 C</b>	<b>Spare filter for Wapo-Ex 8 Plus/ Wapo-Ex 8 III/ Wapo-Ex 8 / Wapo-Ex old up to #40189</b> consisting of 4 filter mats 400 x 325 mm	62,00
<b>611013 C</b>	<b>Spare filter for Wapo-Ex</b> consisting of 10 filter mats 400 x 145 mm	62,00
<b>655001 C</b>	<b>Wax remover, 2 liter</b>	123,00
<b>655002 C</b>	<b>Lime remover, 1000 ml</b>	100,00
<b>!</b>	<b>When you order new filters, please mention the S/N serial number of the unit.</b>	<b>!</b>

### **Compactors**

<b>170993 C</b>	<b>Polyquick 260-T105</b> (combination automat, see page 12)	4.359,00
<b>170976 C</b>	<b>Polyquick 160-T95</b> (combination automat, see page 12)	3.499,00
<b>180994 A</b>	<b>Compactor DEG-1 II</b>	2.390,00
<b>380032 C</b>	<b>O-Ring for DEG / Polyquick</b>	30,55
<b>380063 C</b>	<b>O-Ring for Polyquick 260</b>	77,80

<b><u>Mixing and dosing unit</u></b>		<b>Price in €</b>
<b>172980 A</b>	<b>Sidomix II</b> Dosage and mixing unit for addition-vulcanizing duplicating silicone, 1:1	1.899,00
 <b><u>Accessories for mixing and dosing unit</u></b>		
<b>172640 C</b>	<b>Set material cask, cpl. 2 x 2 liter</b> incl. quick fitting coupling	113,00
<b>172631 C</b>	<b>Stand small for WaSil silicone 2 x 6 kg cask</b>	219,00
<b>172632 C</b>	<b>Stand tall for WaSil silicone 2 x 25 kg cask</b>	293,00
<b>172633 C</b>	<b>Connecting tubes cpl., 2 pcs.</b>	73,50
<b>172634 C</b>	<b>Adapter set for WaSil silicone 2 x 1 kg cask</b>	66,00
<b>172636 C</b>	<b>Adapter set for WaSil silicone 2 x 6 kg cask</b>	73,50
<b>622002 C</b>	<b>Mixing canulas 1:1, 5.4 mm inner Ø, pink, 100 pcs.</b> connectable to dosage automat 1:1 Sidomix II, Sidomix, DA 2000 and similar equipment	81,00
<b>!</b>	<b>For a long durability of the mixing and dosing unit we recommend the use of our Wassermann WaSil duplicating silicone. More on Page 25.</b>	<b>!</b>

	<u>Ejection press</u>	Price in €
<b>128999 A</b>	<b>APR-1</b> for cast rings/ cast flasks (DeguDent/ Heraeus Kulzer), cast ring set 1	1.125,00
	<u>Flask presses</u>	
<b>126999 A</b>	<b>Hydraulic press WW-33 for 1-3 flasks</b>	1.450,00
<b>126998 A</b>	<b>Hydraulic press WW-33 for 2-4 flasks</b>	1.569,00
<b>129997 C</b>	<b>Universal flask clamp for 1 flask</b>	206,00
<b>129996 C</b>	<b>Universal flask clamp for 2 flasks</b>	213,00
<b>129995 C</b>	<b>Universal flask clamp for 3 flasks</b>	218,00
	<u>Accessories for flask presses</u>	
<b>129994 C</b>	<b>Stainless steel clamp clip for 1 flask</b>	155,00
<b>129993 C</b>	<b>Stainless steel clamp clip for 2 flasks</b>	161,00
<b>129992 C</b>	<b>Stainless steel clamp clip for 3 flasks</b>	171,00
<b>129990 C</b>	<b>Universal flask</b>	134,00
<b>629005 C</b>	<b>Spare bottle hydraulic oil, 500 ml</b> for WW-11, WW-22, WW-33 and KÜVOMAT	100,00
<b>629004 C</b>	<b>Hydraulic Oil, Bottle 100 ml</b>	35,50



<u>Surface finish</u>	<b>Price in €</b>
<u>Fine blasting units</u>	
<b>190984 A Cemat-4 II</b> with 4 blasting material containers and sand-drawer, individually adjustable working pressure per blasting chamber	6.290,00
<b>190983 A Cemat-3 II</b> with 3 blasting material containers and sand-drawer, individually adjustable working pressure per blasting chamber	5.650,00
<b>190982 A Cemat-2 II</b> with 2 blasting material containers and sand-drawer, individually adjustable working pressure per blasting chamber	5.029,00
<b>190884 A Gobi Comfort-4</b> with 4 blasting material containers and sand-drawer,	4.699,00
<b>190883 A Gobi Comfort-3</b> with 3 blasting material containers and sand-drawer,	4.099,00
<b>190882 A Gobi Comfort-2</b> with 2 blasting material containers and sand-drawer,	3.450,00
<b>190881 A Gobi Classic</b> with 2 blasting material containers	2.399,00
<b>190885 A Puna II</b> with 2 blasting material containers	1.649,00

**! The fine blasting units are perfectly matched to the combination with our low-pressure suction units. For a long durability of the fine blasting units we recommend the use of Wassermann blasting materials. !**

**Wassermann suction units starting from page 21**

**Wassermann blasting materials starting from page 26.**

**Accessories for fine blasting units**
**Price in €**
**Accessories Cemat II/NT, Gobi Comfort/Classic, Puna/II**

For Cemat OLD please see page 20

<b>190553 C</b>	Spare fine blasting nozzle 0.4 mm (only Cemat-NT)	107,00
<b>190550 C</b>	Spare fine blasting nozzle 0.6 mm	107,00
<b>190551 C</b>	Spare fine blasting nozzle 0.8 mm	107,00
<b>190552 C</b>	Spare fine blasting nozzle 1.2 mm	107,00
<b>190565 C</b>	Spare nozzle (tank) 2.0 mm	18,00
<b>190566 C</b>	Spare nozzle (tank) 1.4 mm	18,00
<b>190494 C</b>	Key for assembly/ disassembly of the tank nozzle	32,50
<b>190175 C</b>	Sleeve for alle versions of Cemat/ Gobi/ Puna, per piece	35,50

**Accessories for Cemat II**

<b>190489 C</b>	Spare hose -yellow-	106,00
<b>190485 C</b>	Spare hose -blue-	106,00
<b>190475 C</b>	Spare hose -red-	106,00
<b>190487 C</b>	Spare hose -green-	106,00
<b>190432 C</b>	Spare protection sheet with special coating	107,00

**Accessories for Gobi Comfort and Gobi Classic**

<b>190739 C</b>	Spare hose -yellow-	100,00
<b>190738 C</b>	Spare hose -blue-	100,00
<b>190729 C</b>	Spare hose -red-	100,00
<b>190728 C</b>	Spare hose -green-	100,00
<b>190432 C</b>	Spare protection sheet with special coating	107,00

**Accessories for Puna II**

<b>190475 C</b>	Spare hose -red-	106,00
<b>190485 C</b>	Spare hose -blue-	106,00
<b>190686 C</b>	Spare protection sheet with special coating	107,00

**Accessories for fine blasting units (previous versions)**      **Price in €**

**Accessories for Cemat-NT (Item no. 190992/ 190990/ 190993)**

Spare Nozzles see page 18

<b>190432 C</b>	Spare protection sheet with special coating from serial number 128538	107,00
<b>190458 C</b>	Conversion kit protection sheet with special coating up to serial number 128537	122,00
<b>190230 C</b>	Spare protection sheet up to serial number 128537	107,00
<b>612003 C</b>	Spare protection foil (plastic), 1 set = 5 pieces, up to serial number 128537	81,00
<b>190239 C</b>	Spare hoses	122,00
<b>190175 C</b>	Sleeve for alle versions of Cemat/ Gobi/ Puna, per piece	35,50

**Accessories for Gobi-2 (Item no. 190896)**

Spare Nozzles see page 18

<b>190432 C</b>	Spare protection sheet with special coating from serial number 128423	107,00
<b>190636 C</b>	Conversion-kit protection sheet with special coating # 124048 - # 127319	387,00
<b>190534 C</b>	Conversion-kit protection sheet for 190899 Gobi-2 up to serial number 124047	339,00
<b>612004 C</b>	Spare protection foil (plastic), 1 set = 5 pieces, up to serial number 127319	81,00
<b>190487 C</b>	Spare hose for green chamber	106,00
<b>190475 C</b>	Spare hose for red chamber	106,00
<b>190485 C</b>	Spare hose for blue chamber	106,00
<b>190489 C</b>	Spare hose for yellow chamber	106,00
<b>190175 C</b>	Sleeve for alle versions of Cemat/ Gobi/ Puna, per piece	35,50

**Accessories for Puna (Item no. 190897)**

Spare Nozzles see page 18

<b>190686 C</b>	Spare protection sheet with special coating	107,00
<b>190475 C</b>	Spare hose for red chamber	106,00
<b>190485 C</b>	Spare hose for blue chamber	106,00
<b>190175 C</b>	Sleeve for alle versions of Cemat/ Gobi/ Puna, per piece	35,50

<b><u>Accessories for Cemat OLD</u></b>		<b>Price in €</b>
<b><u>(Item no. 190996/ 190999/ 190998/ 190997)</u></b>		
<b>190080 C</b>	Spare fine blasting nozzle 0.4 mm	122,00
<b>190081 C</b>	Spare fine blasting nozzle 0.6 mm	122,00
<b>190083 C</b>	Spare fine blasting nozzle 0.8 mm	122,00
<b>190084 C</b>	Spare fine blasting nozzle 1.2 mm	122,00
<b>190111 C</b>	Spare hoses	100,00
<b>612001 C</b>	Spare protection foil, 1 set = 5 pieces	81,00
<b>190050 C</b>	Spare protection sheet (rectangular)	96,50

### **Polishing compact units**

<b>146979 B</b>	<b>Polishing unit WP-Ex 3000 II</b> incl. adjustable motor (750-6000 min <sup>-1</sup> ), dust extraction, LED illumination and LED Spot, rubber polishing troughs, drawer, brush spindle left + right, filters, 2 storage bowls	5.599,00
<b>146978 B</b>	<b>Polishing unit WP-Ex 2000 II</b> incl. 500 W motor, dust extraction, LED illumination, rubber polishing troughs, drawer, brush spindle left + right, filters, 2 storage bowls	4.559,00
<b>146971 B</b>	<b>Polishing unit WP-Ex 10 II</b> incl. 500 W motor, dust extraction, LED illumination, rubber polishing trough, brush spindle right , filter, 1 storage bowl	2.990,00

### **Polishing motor**

<b>155999 B</b>	<b>WSM-2 (1500/ 3000 min<sup>-1</sup>)</b> incl. brush spindle left + right	1.075,00
-----------------	--	----------

### **Polishing troughs**

<b>143123 C</b>	<b>Polishing trough without extractor nozzle (2017)</b> incl. adjustable protection sheet and removable rubber trough	253,00
<b>143115 C</b>	<b>Polishing trough with extractor nozzle (2017)</b> (Ø75 mm), incl. cap incl. adjustable protection sheet and removable rubber trough	287,00
<b>143131 C</b>	<b>Polishing trough with extractor nozzle (2017)</b> (Ø40 mm), incl. cap incl. adjustable protection sheet and removable rubber trough	287,00
<b>143550 C</b>	<b>Extraction traverse II, with switch to guide the air</b>	568,00

output nozzle Ø75 mm

<u>Accessories for polishing compact units/ polishing motors</u>		Price in €
<b>152025 C</b>	Brush spindle 92 mm for left side	54,00
<b>152022 C</b>	Brush spindle 92 mm for right side	54,00
<b>152057 C</b>	Stone holder for left side	54,00
<b>152056 C</b>	Stone holder for right side	54,00
<b>152044 C</b>	Brush spindle 125 mm for left side	81,00
<b>152043 C</b>	Brush spindle 125 mm for right side	81,00
<b>152049 C</b>	Chuck clamping from 1-6 mm	212,00
<b>155123 C</b>	Wrench for pull nut	65,50
<b>146390 C</b>	Daylight LED Spot	286,00
<b>146585 C</b>	<b>RepairSet Glass Panel from old to new</b>	458,00

Filter for polishing compact units

<b>611005 C</b>	<b>Spare filter bag for WP-Ex 3000 II/ 2000 II/ 10 II</b> dust class M, for micro fine dust, 1 set = 2 pieces	95,50
<b>611011 C</b>	<b>Spare filter bag WP-Ex 3000/ 2000/ 10</b> dust class M, for micro fine dust, 1 set = 2 pieces	81,00
<b>611021 C</b>	<b>Spare filter bag for WP-Ex 80</b> dust class M, for micro fine dust, 1 set = 2 pieces	81,00

**!** Please sent us the serial number of your unit in case of ordering filter material. The serial number allows us to find quickly without any doubt the correct sparteparts of your unit.  
**!** The identification plate is located on the back or right side of your unit.

Automatic stand

<b>146170 D</b>	<b>Automatic stand for WP-Ex 3000 II/ 2000 II</b>	2.990,00
	vertically adjustable from 74 – 109 cm Dimension (table top): D 560 x W 675 mm	

<u>Suction technology</u>		Price in €
<b>140987 A</b>	<b>Suction unit SG-10 activated carbon version</b> incl. filter and module with activated carbon filter Available on request, can be retrofitted: <b>140408 C</b> SG-10 External speed regulator for fan control	2.169,00
<b>140988 A</b>	<b>Suction unit SG-10 standard version</b> incl. filter, without accessories Available on request, can be retrofitted: <b>140408 C</b> SG-10 External speed regulator for fan control	1.385,00
<b>140997 A</b>	<b>Single-place-dust extraction SG-1/2 D</b> equipped with automatic switch-off in case of exhausted filter capacity, incl. filter, without accessories	1.835,00
<b>140970 C</b>	<b>Suction unit SG-Economy <span style="color: red;">New!</span></b> incl. 3 stage filter system, without accessories	1.350,00
<b>140995 A</b>	<b>Single place dust extraction SG-1/1</b> incl. filter, without accessories	1.099,00
<b>143100 C</b>	<b>Automatic starting device AS-100</b> for automatic connection of extractor systems	422,00

	<b>Price in €</b>
<b>145997 A Dust extraction drawer LSG-02 DE</b> complete with universal mounting frame and automatic switch-off in case of exhausted filter capacity, incl. filter , additional equipped with auto-start-system	2.330,00
<b>145998 A Dust extraction drawer LSG-02 D</b> complete with universal mounting frame and automatic switch-off in case of exhausted filter capacity, incl. filter	2.150,00
<b>145240 C Activated charcoal filters for monomer for LSG-02 DE/ LSG-02 D</b> for elimination of dangerous vapors	668,00
<b>145200 C Universal extractor hood for LSG-02 D and DE</b> with exchangeable protection sheet	400,00
<b>145005 C Additional drawer for LSG-02 D and DE</b> can also be retrofitted	216,00
<b>149999 A Mobile dust extraction Timobil</b> cpl. with handpiece motor 1.000 – 40.000 min <sup>-1</sup> , 160 Watt, infinitely adjustable, standard height 900 mm	4.559,00
<b>149130 C Extractor hood for Timobil</b> Closed, can also be retrofitted	391,00
<b>149140 C Extractor hood illuminated for Timobil</b>	599,00
 <b><u>Suction Box</u></b>	
<b>143390 C Suction Box Compact III</b> with daylight LED illumination and 1 extractor nozzle Ø 75 mm exterior (optional mounting in the front or in the rear)	789,00
 <b><u>Accessories for Suction Box Compact III + previous versions</u></b>	
<b>143303 C</b> Set of arm rests	81,00
<b>143305 C</b> Set of cover plates	210,00

190686 C Spare protection sheet with special coating 107,00

**Accessories for suction units**

**Price in €**

**Telescopic channels + corresponding suction parts**

143070 C	Telescopic channel 570 mm	144,00
143090 C	Telescopic channel 455 mm	158,00
143074 C	Telescopic channel 480 mm	153,00
143066 C	Suction part, transverse	114,00
143080 C	Suction part, straight	125,00
145245 C	Activated charcoal filters for monomer - suction part, straight for elimination of dangerous vapors	668,00
143108 C	Cap for suction part transverse	48,50
143200 C	Universal hood for suction part straight	400,00
143220 C	Suction part with wooden rest	338,00
143221 C	Suction part with wooden rest/plastic screen	401,00
143231 C	Suction part with wooden rest/glass screen	447,00

**Special channel + appropriate suction parts**

(e.g. for installation in Freuding tables)

143093 C	Special channel with hose coupling	174,00
143041 C	Extractor nozzle	144,00

**! Channel and one suction part are always both required!**

**!**

**Hoses**

140057 C	PVC-hose 75 Ø - 1.20 m length, incl. 2 hose clips	107,00
320011 C	PVC-hose 75 Ø - 1.80 m length	152,00
320010 C	PVC-hose 75 Ø - 2.40 m length	187,00
320012 C	PVC-hose 75 Ø - 3.60 m length	262,00
322003 C	Hose clip 75 Ø	5,10
327023 C	PVC-hose 51 Ø inner - per meter	58,50

**Special accessories**

140056 C	Throttle valve 75 Ø	162,00
----------	---------------------	--------

**! Reduction- and transition pieces of all measurements available.**

**!**

**Would you please be so kind as to mention inner- and outer diameter in case of supplementary order.**



<u>Filters for suction units</u>		Price in €
<b><u>Filters &amp; accessories for suction unit SG-10</u></b>		
<b>611010 C</b>	Spare filter bag yellow dust class M, for micro fine dust, 1 set = 4 pieces	191,00
<b>611018 C</b>	Spare filter bag white, for plaster class 1 and 2 only (!) regenerative, 1 set = 4 pieces	191,00
<b>610073 C</b>	*Filter bag white (plaster class 1/ class 2), big	75,50
<b>610074 C</b>	*Filter bag yellow (micro fine dust), big <i>* to be used only in combination with item no. 140410</i>	75,50
<b>140410 C</b>	Filter extension module SG-10	258,00
<b>140435 C</b>	Filter cassette, cpl. SG-10/ SG-20	88,50
<b>140424 C</b>	Module activated carbon filter SG-10	513,00
<b>140436 C</b>	Activated carbon filter (AKM 6-ply) SG-10	150,00
<b><u>SG-1/1 + SG-2/2 + Timobil</u></b>		
<b>611010 C</b>	Spare filter bag yellow, dust class M, for micro fine dust, 1 set = 4 pieces	191,00
<b>611018 C</b>	Spare filter bag white, (class 1/class 2) for plaster, regenerable, 1 set = 4 pieces	229,00
<b><u>SG-1/2 D+ LSG-02 DE/ LSG-02 D</u></b>		
<b>611020 C</b>	Spare filter bag yellow, dust class M, for micro fine dust, 1 set = 4 pieces	229,00
<b>611017 C</b>	Spare filter bag white, (class 1/class 2) for plaster, regenerable, 1 set = 3 pieces	173,00
<b>610056 C</b>	Spare filter monomer for LSG-02	69,00
<b>611022 C</b>	<b><u>Filter-Kit SG-Economy</u></b> 1 x Spare-main filter / cartridge filter SG-Economy 5 x Spare-filter bag SG-Economy 2 x Spare-HEPA exhaust filter SG-Economy	182,00
<b><u>Filters &amp; accessories for suction unit SG-20</u></b>		
<b>610073 C</b>	Filter bag white (plaster class 1/ class 2)	75,50
<b>610074 C</b>	Filter bag yellow (micro fine dust)	75,50
<b>140551 C</b>	Activated carbon filter, cpl. SG-20	145,00
<b>140435 C</b>	Filter cassette, cpl. SG-10/ SG-20	88,50
<b>610019 C</b>	Motor protective filter SG-20	11,50
<b><u>WTS-1R</u></b>		
<b>611019 C</b>	filter bag FB 25, 1 Set = 3 pieces	107,00
<b>!</b>	<b>Filters for previous versions of our suction units on demand.</b>	<b>!</b>

Consumable materials

Price in €

Duplicating silicone
**WaSil soft**

Addition-vulcanizing duplicating silicone, 1:1  
 hardness Shore A: 16-18  
 component A white (catalyst)  
 component B yellow (base)

<b>656060 C</b>	<b>WaSil soft 2 x 1 kg component A/ B</b>	110,00
<b>656061 C</b>	<b>WaSil soft 2 x 6 kg component A/ B</b>	565,00
<b>656062 C</b>	<b>WaSil soft 2 x 25 kg component A/ B</b>	2.290,00

**WaSil speed**

Addition-vulcanizing duplicating silicone, 1:1  
 high density, hardness Shore A: > 24  
 component A white (catalyst)  
 component B turquoise (base)

<b>656070 C</b>	<b>WaSil speed 2 x 1 kg component A/ B</b>	125,00
<b>656071 C</b>	<b>WaSil speed 2 x 6 kg component A/ B</b>	696,00
<b>656072 C</b>	<b>WaSil speed 2 x 25 kg component A/ B</b>	2.710,00

Mixing canulas, 1:1

5.4 mm inner Ø

<b>622002 C</b>	<b>mixing canulas, pink, 100 pcs.</b> connectable to dosage automat 1:1 Sidomix II, Sidomix, DA 2000 and similar equipment	81,00
-----------------	--	-------

Price in €

**Blasting materials****Aluminium oxyde**

For blasting of all dental alloys. The purity, sharp-edged grains and extreme hardness guarantee very efficient blasting.

Silicosis-free, made in Germany

<b>656010 C</b>	<b>Aluminium oxyde 50 µm – 10 kg</b>	90,50
<b>656015 C</b>	<b>Aluminium oxyde 50 µm – 20 kg</b>	187,00
<b>656013 C</b>	<b>Aluminium oxyde 110 µm - 10 kg</b>	90,50
<b>656016 C</b>	<b>Aluminium oxyde 110 µm - 20 kg</b>	187,00
<b>656014 C</b>	<b>Aluminium oxyde 250 µm - 10 kg</b>	90,50
<b>656018 C</b>	<b>Aluminium oxyde 250 µm – 20 kg</b>	187,00

**Glass beads**

For gentle cleaning and compacting of sensitive surfaces of precious and non-precious alloys. For shining and smoothing out of crowns, bridges and CrCo techniques.

Silicosis-free, made in Germany

<b>656020 C</b>	<b>Glass beads 50 µm – 8 kg</b>	70,50
<b>656021 C</b>	<b>Glass beads 50 µm – 15 kg</b>	99,50
<b>656019 C</b>	<b>Glass beads 125 µm – 8 kg</b>	70,50
<b>656022 C</b>	<b>Glass beads 125 µm – 15 kg</b>	99,50



Wassermann Dental-Maschinen GmbH  
Rudorffweg 15–17 · 21031 Hamburg · Germany  
Phone: +49 (0)40 730 926-0  
[export@wassermann-dental.com](mailto:export@wassermann-dental.com) · [www.wassermann.hamburg](http://www.wassermann.hamburg)



[/WassermannDentalmaschinen](https://www.facebook.com/WassermannDentalmaschinen)