

“developed  
for more  
**flexibility**”

- **efficient & powerful**

**Our Proven Bestseller:**

Wassermann's best-selling trimmer  
- now even stronger, even more  
efficient, even more flexible.

**OUR PROVEN RECIPE FOR SUCCESS:**

We combine technology, experience,  
competence and team spirit.



**Always available**

Would you like to request a quote or  
place an order? Or do you have  
questions about your order, delivery  
or warranty?

We look forward to your call and  
will be happy to help you.

**Dry Trimmer  
HSS-TR**



**Service Hotline**

+49 (0)40 730 926-0



[www.wassermann.hamburg](http://www.wassermann.hamburg)



[fb.com/WassermannDentalmaschinen](https://fb.com/WassermannDentalmaschinen)

SCAN ME 



# Dry Trimmer HSS-TR

**NEW** Powerful  
in a new look

  
QUALITY MADE IN GERMANY SINCE 1927

Its unique, central connection for the suction system guarantees the best suction performance. Due to its powerful, high-torque but also energy-efficient motor highest grinding speeds are achieved.



## The most functions and features at a glance:

- Versatile: Various Wassermanns grinding wheels can be used
- Very high-torque and maintenance-free motor (at max. 3000 rpm)
- Stable and solid stainless steel and noble aluminum
- Central connection for suction system
- Powerful disc suction and internal suction
- Safe and easy to open, no unnecessary seals
- Large work area of 10,5 x 13,5 cm
- Grinding jigs KMS are also available in Swiss and EOS type.

Notes:

### Technical Data:

Output: 1.000 W  
Power consumption: 4,6 A

Weight without grinding disc: 19,8 kg  
W x H x D 330 x 335 x 410 mm

Speed: 3000 U/min  
Voltage: 220–240 V

### Support Service:

Wassermann Dental-Maschinen GmbH  
Rudorffweg 15 - 17  
21031 Hamburg

Tel: +49 (0)40 730 926-0  
Fax: +49 (0)40 730 37 24  
www.wassermann-dental.com  
sales@wassermann-dental.com